# What you should know about Juzo Compression Garments.





For a healthier, more active life

...Keep Moving.









For over 100 years, it's been our pleasure to provide compression therapy garments that help manage symptoms without inhibiting a person's lifestyle.

### Frequently Asked Questions

#### When to wear garment?

First thing in the morning and remove at the end of the day.

# When should I purchase a new compression garment?

Every 6 months or as therapy progresses and new symptoms develop.

#### Can I use lotions & ointments? What about Latex allergies?

Juzo® compression garments are LATEX FREE using high quality Lycra® covered with textile fibers that offer wearing comfort, air permeability, and are largely resistance to sun, ointments, lotions and perspiration.

# Where can I find material content, size and style?

A tag sewn inside provides details of the garment.

#### What's inside a Juzo compression package?

Care & application brochure highlight use and related information. Open toe stockings include a donning device "open toe slippie" and instructions.

## What does graduated compression mean?

Compression begins at the ankle and decreases up the leg. For arm sleeves, the compression begins at the wrist and decreases moving up the arm.

#### Care and Maintenance:

# How to wash your garment?

Washing helps restore the garment back to form.

Hand or machine wash warm daily.

Use the lowest permanent press setting.

Garment bag is recommended.

Separate dark colors from other clothing.

Use Juzo® washing detergent or another mild soap.

DO NOT USE: bleach, chlorine, fabric softeners, or other laundry additives.

Always use Juzo latex free donning gloves when applying compression garments. Rough finger nails, jagged rings etc. can damage fabric.

Watch a video demonstration on how to take care of your compression garment & donning/doffing a compression garment at the Juzo websites below:

www.juzousa.com

Patient Therapy Line: 888-255-1300 • support@juzousa.com