



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

## FarrowWrap® Prescription Form

### Prescribing Information



	Left	Right
Thigh	<input type="checkbox"/>	<input type="checkbox"/>
Leg	<input type="checkbox"/>	<input type="checkbox"/>
Foot	<input type="checkbox"/>	<input type="checkbox"/>
Toes	<input type="checkbox"/>	<input type="checkbox"/>



Mild to Moderate Edema

Moderate to Severe Edema\*

\*Be cautious when prescribing "Moderate to Severe Edema" for patients with PAD.

If mixed severity, please specify in special instructions below.

	Left	Right
Arm	<input type="checkbox"/>	<input type="checkbox"/>
Hand	<input type="checkbox"/>	<input type="checkbox"/>
Fingers	<input type="checkbox"/>	<input type="checkbox"/>




Special Instructions: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

Related Diagnoses: \_\_\_\_\_

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

# FarrowHybrid® Knee High OTS and Custom Order Form

Note: All measurements in cm



Please note that “\_” in the SKU below (i.e. FW\_-O-LR) stands for either “CL” for Classic, “LT” for LITE, “ST” for STRONG or “BA” for BASIC. The SKU's may then be cross-referenced with our price sheet to obtain the garment price.

Please place a numeric value in the white box to indicate items needed.		Off-The-Shelf Legpiece Size Chart									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(C) Widest Calf		36 – 43 cm		42 – 50 cm		48 – 58 cm		53 – 63 cm		58 – 68 cm	
(B) Ankle		21 – 25 cm		25 – 30 cm		30 – 36 cm		36 – 42 cm		42 – 50 cm	
(A-D) Reg	FW_-O-LR	33 – 37 cm		35 – 39 cm		37 – 41 cm		39 – 43 cm		39 – 43 cm	
(A-D) Tall	FW_-O-LT	38 – 41 cm		40 – 43 cm		42 – 45 cm		44 – 47 cm		44 – 47 cm	

Please place a numeric value in the white box to indicate items needed.		Off-The-Shelf Footpiece Size Chart									
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty
(A <sup>1</sup> ) Mid-Foot		22 – 24 cm		25 – 27 cm		28 – 30 cm		31 – 34 cm		35 – 40 cm	
(X) Reg	FW_-O-FR	16 – 17 cm		17 – 18 cm		19 – 20 cm		20 – 21 cm		22 – 23 cm	
(X) Long	FW_-O-FL	18 – 19 cm		19 – 20 cm		21 – 22 cm		22 – 23 cm		24 – 25 cm	

Options	
<b>FABRIC</b>	✓
Classic	
LITE	
STRONG	
BASIC*	
Liners	
Small Silver Liner < 60 cm	
Large Silver Liner < 70 cm	
Terry Cloth Liner 30 – 120 cm	

Custom Components	Legpiece	FW_-C-L		Footpiece	FW_-C-F	
Quantity	Right:		Left:	Right:		Left:

\*BASIC is available as Off-The-Shelf only.

When a footpiece & legpiece are purchased together, 1 pair of liners are included free. Additional liners will cost extra. See liner form & pricing for details. *BASIC foot and leg piece combinations will only be issued one (1) line; not a pair.*

# Directions for Knee High Garments

Please note that OTS sizes can accommodate at least a 15 to 20% reduction in edema, but only a 5% increase.

Fill in the Patient / Billing & Shipping Information areas **completely**. Next, choose the FarrowWrap® version; Classic, LITE, STRONG or BASIC. Classic, STRONG and BASIC are generally reserved for more severe edema cases. LITE is more appropriate for mild to moderate edema. On the measurement page, part of the SKU number (i.e. FW\_-O-LR) is represented by a “\_” which represents either a “CL” for Classic, “LT” for LITE, “ST” for STRONG or “BA” for BASIC. (**NOTE: BASIC is not available as a Custom**) The appropriate SKU number can be cross-referenced with the Pricing Schedule to calculate the price.

Only circumferences at points **A**, **B**, & **C** and length measurements **X** & **A-D** are necessary to determine if the patient can fit into an OTS (Off-The-Shelf (ready-made) garment. If the patient will not fit into an OTS garment, then the additional measurements for that particular component are required. Length measurements are measured along the posterior contour of the leg following contours around any lobes. If any lobes are present, drawings and/or pictures (with the patient’s permission) will help. Measure only portions of the limb to be fitted with a garment. Follow these instructions on how to properly measure. ONLY record ACTUAL measurements.

## FarrowWrap® Legpiece

1. Obtain the following 3 measurements:
  - a. Widest calf circumference at **C**.
  - b. Least ankle circumference **B**, just above malleoli.
  - c. Length of leg following posterior contour from 2 finger widths below posterior knee crease to the ground, **A-D**.
2. Match the OTS widest calf and ankle circumferences to a respective size. If the mid-calf AND ankle circumferences do not both match a single size, or if the patient has a large lobule, then a custom legpiece is highly recommended and we will need all associated measurements. (\*BASIC is not available as a custom)
3. If the circumferences do match an OTS size, next compare the **A-D** length to the OTS **Regular** and **Tall** measurement ranges. If the **A-D** falls within one of these ranges, then in the space next to the appropriate OTS length enter the number of legpieces needed in that given size. If the measured length is not within the length range of a Regular or Tall, a custom legpiece is highly recommended and all associated measurements would be required. Enter the number of custom legpieces needed at the bottom of the measurement form in the **Custom Components** section.

## Note

The leg measurements used by this form allow for at least a 20% additional reduction in leg and ankle circumference (edema reduction), but can only accommodate about a 5% increase. There is room for greater variance in calf circumferences than for ankle circumferences. Patients may fit into an OTS garment even though this form may indicate otherwise, but we cannot guarantee it.

## FarrowWrap® Footpiece

1. Measure the affected foot’s mid-foot circumference, **A**, and length, **X**. The length is measured from the proximal (posterior) border of either the 1st or 5th metatarsal head (bunion), which ever is shorter, to the posterior most aspect of the heel. This measurement should be made across the floor - not along the side of the patient’s foot.
2. Next, match the **A** circumference to an OTS footpiece size. If a match cannot be made, a custom footpiece will be required.
3. Next, compare the **X** length to the OTS **Regular** and **Long** length measurement ranges. If the **X** length falls within one of these OTS lengths, then in the space next to the appropriate OTS length enter the number of footpieces you would like in that size. If the measured length is not within the length range of a Regular or Long, a custom footpiece is highly recommended. Enter the number of custom footpieces needed at the bottom of the measurement form in the **Custom Components** section.

Note: The **X** measurement defaults to extend to the metatarsal heads for safety reasons. If the footpiece needs to extend further out, then simply measure the straight distance to the desired point. The footpiece is not meant to be a "closed toe" item. **Only do this on patients with adequate circulation and no peripheral neuropathy; do not use on symptomatic diabetics.**

\*As a replacement option in mild to moderate edema, a FarrowHybrid® AD Foot Compression may be used in place of a footpiece.

## Liner Options

Unless specified, a pair of **Small Farrow Silver Liners (FLSF-AD2)** will be sent **Free** with each Extra Small and Small footpiece / legpiece combination and a pair of **Large Farrow Silver Liners (FLSL-AD4)** will be sent with each Medium, Large and Extra Large footpiece / legpiece combination. **FarrowWrap® BASIC** garments will follow the same format but will only be issued one (1) liner, not a pair. The FLSF-AD2 can accommodate leg circumferences up to 60 cm and the FLSL-AD4 can accommodate leg circumferences up to 70 cm. If the leg circumference is greater than 70 cm, a pair of **Terry Cloth AD Liners** will be issued instead. Terry Cloth Liners come in 3 sizes; S < 40 cm, M 40-70 cm, & L 70-125 cm. They may be **substituted** for silver liners in the case of silver incompatibilities without charge. Indicate on the order form with a check mark which liner is required. To order additional liners use the liner order form or indicate the amount needed on the "**Notes**" line at the bottom of the page. Additional liners will be billed. Please note if **ordering a FarrowFoam™ Liner for use under an OTS garment, that you must also add 10 cm to leg circumferences and 5 cm to foot circumferences when determining the correct sized FarrowWrap® leg or foot piece.**

## Other Options


Some patients have very triangular shaped legs and even the Classic garment may have difficulty staying up. Additional velcro-like spines, called **Hook Stays**, can be easily attached to the garment to help hold it in place. The stays are available in small (16 cm length), medium (21 cm length), and large (27 cm length). For custom garments, additional velcro tabs may be ordered to help prevent gapping between bands on unusually shaped limbs.

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

# FarrowWrap® Thigh High OTS and Custom Order Form



Please place a numeric value in the white box to indicate items needed.											<b>Off-The-Shelf Thighpiece Size Chart</b> (Only comes with wrap-around kneepiece)										
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty										
(G) Groin		60-70 cm		65-75 cm		70-80 cm		75-90 cm		85-100 cm											
(E) Above Knee		47-52 cm		53-59 cm		60-66 cm		67-73 cm		74-80 cm											
(E-G) Short	FW_-O-TS	16-20 cm		16-20 cm		16-20 cm		16-20 cm		16-20 cm											
(E-G) Reg	FW_-O-TR	21-26 cm		21-26 cm		21-26 cm		21-26 cm		21-26 cm											
(E-G) Tall	FW_-O-TT	27-31 cm		27-31 cm		27-31 cm		27-31 cm		27-31 cm											
Please place a numeric value in the white box to indicate items needed.											<b>Off-The-Shelf Legpiece Size Chart</b>										
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty										
(C) Widest Calf		36-43 cm		42-50 cm		48-58 cm		53-63 cm		58-68 cm											
(B) Ankle		21-25 cm		25-30 cm		30-36 cm		36-42 cm		42-50 cm											
(A-D) Reg	FW_-O-LR	33-37 cm		35-39 cm		37-41 cm		39-43 cm		39-43 cm											
(A-D) Tall	FW_-O-LT	38-41 cm		40-43 cm		42-45 cm		44-47 cm		44-47 cm											
Please place a numeric value in the white box to indicate items needed.											<b>Off-The-Shelf Footpiece Size Chart</b>										
		XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty	XLarge	Qty										
(A <sup>1</sup> ) Mid-Foot		22-24 cm		25-27 cm		28-30 cm		31-34 cm		35-40 cm											
(X) Reg	FW_-O-FR	16-17 cm		17-18 cm		19-20 cm		20-21 cm		22-23 cm											
(X) Long	FW_-O-FL	18-19 cm		19-20 cm		21-22 cm		22-23 cm		24-25 cm											

Options	
<b>FABRIC</b>	<input checked="" type="checkbox"/>
Classic	<input type="checkbox"/>
LITE	<input type="checkbox"/>
STRONG	<input type="checkbox"/>
When a footpiece, legpiece and thighpiece are purchased together, 1 pair of <b>Terry Cloth</b> liners are included automatically free. Other liners may cost extra. See liner form & pricing for details.	

Custom Component	Thighpiece	<input type="checkbox"/> Wrap-around Standard	Legpiece	FW_-C-L	Footpiece	FW_-C-F	Extend footpiece length to base of toes
Quantity	Right:	Left:	Right:	Left:	Right:	Left:	L R BL:

Other: \_\_\_\_\_

# Directions for Thigh High Garments

Please note that OTS sizes can accommodate about a 20% reduction in edema, but only a 5% increase.

Fill in the Patient / Billing & Shipping Information areas **completely**. Next, choose the FarrowWrap® version; Classic, LITE or STRONG. Classic and STRONG are generally reserved for more severe edema cases. LITE is more appropriate for mild to moderate edema. On the measurement page, part of the sku number is represented by a “\_” which represents either a “CL” for Classic, “LT” for LITE or “ST” for STRONG. The appropriate SKU number can be crossreferenced with the Pricing Schedule to calculate the price.

Only circumferences at points **A**, **B**, **C**, **E**, & **G** and length measurements **X**, **A-D**, & **E-G** are necessary to determine if the patient can fit into an **OTS** (Off-The-Shelf or ready-made) garment. Length measurements for the legpiece are measured along the posterior contour of the leg and length measurements for the thighpiece are measured medially and following contours around any lobes. If any lobes are present, drawings and/or pictures (with the patient’s permission) will help. Follow these instructions on how to properly measure.

## FarrowWrap® Thighpiece

1. Obtain the following 3 measurements:
  - a. Upper thigh circumference (at **G**) just below the groin.
  - b. Distal thigh circumference (at **E**) right above the knee.
  - c. Length of thigh from just above the knee joint to the groin. Measure medially. **E-G**.
2. Now, match the upper and distal thigh circumferences to a respective size. If the upper thigh AND distal thigh circumferences do not both match a single size, or if the patient has a large lobule, then a custom thighpiece is necessary along with all measurements.
3. If the circumferences do match an OTS size, then compare **E-G** to the OTS Short, Regular and Tall length measurements. If the **E-G** length falls within the specified ranges, then in the space next to the appropriate OTS length enter the number of thighpieces you would like in that given size. If the measured length is not within the ranges, a custom thighpiece will be necessary. Enter the number of custom thighpieces required at the bottom of the page in the **Custom Component** section.
4. All OTS thighpieces come with a wrap-around kneepiece. Custom thighpieces come with a standard kneepiece or a wrap-around kneepiece. See our catalog for descriptions. If choosing a custom thighpiece, enter the number of pieces needed for the RIGHT and/or LEFT along with the kneepiece desired, Wrap-Around or Standard.
5. We recommend biker’s shorts or GarmentGrip™ to help hold thigh pieces in place on problematic patients.

## FarrowWrap® Legpiece and Footpiece

1. See FarrowWrap® A-D (Knee High) Order Form Instructions

## Note

The measurements used by this form allow for at least a 20% additional reduction in leg and thigh circumferences (edema reduction), but can only accommodate about a 5% increase. There is room for greater variance in calf circumferences than for ankle circumferences. Patients may fit into an OTS garment even though this form may indicate otherwise, but we cannot guarantee it.

## Liner Options

All footpiece/legpiece/thighpiece combinations will be issued a free pair of Terry Cloth AG Liners automatically. To order additional liners, please use the “Liners Order From” section. We recommend at least 2 pairs of liners per extremity. Terry Cloth Liners come in 3 circumferences, S <40cm, M 40-70 cm, & L 70-120 cm.

## Other Options


Hook Stays can be used to help prevent FarrowWrap® legpieces from sliding on highly triangular legs and can be easily attached to the garment. Stays are available in small (16 cm length), medium (21 cm length), and large (27 cm length).

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowWrap® Off-The-Shelf Arm



FarrowWrap® OTS Arm			Small	Medium	Large
(G) axilla circumference			22 – 31 cm	29 – 39 cm	32 – 45 cm
(E) elbow crease circumference			20 – 27 cm	25 – 34 cm	30 – 40 cm
(C) wrist crease circumference			14 – 18 cm	16 – 21 cm	19 – 25 cm
(C-G) Short	<b>40 – 43 cm</b>	Qty			
(C-G) Regular	<b>44 – 47 cm</b>	Qty			
(C-G) Long	<b>48 – 51 cm</b>	Qty			

Options	
<b>SIDE</b>	✓
Left	
Right	
<b>Color</b>	✓
Tan	
Black	

**Left or Right must be selected in the “Options” box.  
Liners will not be issued with this garment.**

*If no color is selected, Tan is the default option.*

### Directions for the Off-The-Shelf Arm

For length, measure from the outside flexor surface of the wrist to the axilla with a slight bend in the arm. Measure circumference at the wrist, elbow crease and axilla. Compare measurements to the size chart to determine correct size. All circumference measurements must correspond to one off-the-shelf size in order for the garment to fit. This garment can be donned and doffed with one hand. **Liners will not be issued with this garment. The OTS Arm is available in LITE fabric only.**


Other: \_\_\_\_\_

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowHybrid® Custom Arm Order Form



Options	
<b>FABRIC</b>	✓
Classic	
LITE	
STRONG	
Custom arm wraps come <b>standard</b> with 1 pair of Terry Cloth Liners or Silver Arm Liners free. Additional liners cost extra.	

Custom Garment	Armpiece		FW_-C-A
Quantity	Right:		Left:

**This garment requires 2 hands to don and 1 to doff.**

### Directions for the Upper Extremity

Choose the FarrowWrap® version: Classic, LITE or STRONG. Classic and STRONG are generally reserved for more severe edema cases. LITE is more appropriate for mild to moderate edema. On the armpiece table, the “\_” is part of the sku number & stands for either a “CL” for Classic, “LT” for LITE or “ST” for STRONG. This SKU number can be crossreferenced with the Pricing Schedule to calculate the price.

For lengths, measure from the outside flexor surface of the wrist up the arm to the axilla. A pair of Terry Cloth liners is shipped free with each order. To order additional liners, enter the quantity of pairs on the line below.


Other: \_\_\_\_\_

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowHybrid® Foot Compression Order Form



**FarrowHybrid® AD Foot Compression** is part compression stocking (over the foot and ankle) and part sock (ankle to knee) and may be used in place of our footpieces with our legpiece. Easy to apply and remove, they are comfortable and readily fit in patients' shoes.

FarrowHybrid® AD		SKU	Small	Medium	Large	XLarge
Male Shoe Size			3 – 6	7 – 9	10 – 12	13 – 15
Female Shoe Size			4 – 7	8 – 10	11 – 13	14 – 16
<b>FarrowHybrid® ADI</b>	FLFH-ADI	Qty		Qty		Qty
<b>FarrowHybrid® ADI Wide</b>	FLFH-ADIW	Qty		Qty		Qty
<b>FarrowHybrid® ADIIS</b>	FLFH-ADIIS	Qty		Qty		Qty
<b>FarrowHybrid® ADIIS Wide</b>	FLFH-ADIISW	Qty		Qty		Qty

Color Options	
Taupe	<input type="checkbox"/>
Black	<input type="checkbox"/>

*Black available only in FarrowHybrid ADI Wide. If no color is selected, Taupe is the default option.*


FarrowHybrid® Foot Compression Size Selector		
	20 – 30mmHg	30 – 40mmHg
A¹ Circumference 20 – 27cm (standard)	FLFH-ADI	FLFH-ADIIS
A¹ Circumference 28 – 35cm (wide)	FLFH-ADIW	FLFH-ADIISW

FarrowHybrids cover a wide range of circumferences in standard and wide sizes. The maximum leg circumference must be less than 60cm. FarrowHybrid® ADI and ADIW provide Class I (20-30mmHg) compression and FarrowHybrid® ADIIS and ADIISW provide Class II (30-40mmHg) compression. The FarrowHybrid® ADIIS foot compression CONTAINS SILVER and should not be used on patients allergic to silver.

**FarrowHybrid® AD Foot Compression**

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

## Farrow Microfine and FarrowWrap® Gauntlet Order Form



Farrow Microfine Gloves and Toe Caps provide Class I compression and are ambidextrous. The fingers and toes of the garments can be trimmed using scissors and won't fray or roll.

FarrowWrap® LITE Hand Gauntlet is ambidextrous and can be worn with the supplied FarrowFoam™ squares (one 4mm & one 8mm). They are suitable for mild to severe edemas. The gauntlet can be worn by itself or over a Microfine glove, for extra dorsal hand compression.

Microfine Class I Glove		XSmall		Small		Medium		Large	
(A) Base of Fingers		16 – 18 cm		18 – 20 cm		20 – 22 cm		22 – 24 cm	
(B) Fold of Thumb		17 – 19 cm		19 – 21 cm		21 – 23 cm		23 – 26 cm	
(C) Wrist		15 – 18 cm		17 – 20 cm		19 – 22 cm		22 – 24 cm	
<b>Microfine Glove - ½ Finger Length - Black</b>	FSMF-T-HHB	Qty		Qty		Qty		Qty	
<b>Microfine Glove - ½ Finger Length - Copper</b>	FSMF-T-HHC	Qty		Qty		Qty		Qty	
<b>Microfine Glove - ½ Finger Length - Sand</b>	FSMF-T-HHS	Qty		Qty		Qty		Qty	
<b>Microfine Glove - Full Finger Length - Black</b>	FSMF-T-HFB	Qty		Qty		Qty		Qty	
<b>Microfine Glove - Full Finger Length - Copper</b>	FSMF-T-HFC	Qty		Qty		Qty		Qty	
<b>Microfine Glove - Full Finger Length - Sand</b>	FSMF-T-HFS	Qty		Qty		Qty		Qty	
FarrowWrap® Hand Gauntlet		XSmall		Small		Medium		Large	
<b>LITE Ambidextrous Hand Gauntlet</b>	FWLT-O-HGAT	Qty		Qty		Qty		Qty	
<b>Gauntlet Foam &amp; Velcro Replacement Pack</b>	FWLT-O-HGRP	Qty				Qty			

Microfine Class I Toe Cap		XSmall		Small		Medium		Large	
(A) Base of Toes		21 – 23 cm		23 – 25 cm		25 – 27 cm		27 – 29 cm	
<b>Microfine Toe Cap – Black</b>	FSMF-T-MB	Qty		Qty		Qty		Qty	
<b>Microfine Toe Cap – Copper</b>	FSMF-T-MC	Qty		Qty		Qty		Qty	
<b>Microfine Toe Cap – Sand</b>	FSMF-T-MS	Qty		Qty		Qty		Qty	

## Directions for Glove, Gauntlet, and Toe Cap

Fill in the Patient / Billing & Shipping Information areas **completely**. Measure hand and/or foot circumferences at specified areas and record on form. Cross reference measurements with measurements on the appropriate size chart to determine the correct size garment. In the blank white box indicate the number of garments desired. Microfine and FarrowWrap® Hand Gauntlets are sold as each.

### Microfine

#### Microfine Glove Measuring Instructions

Measure circumference at base of fingers, fold of thumb, and wrist. Match measurements to appropriate sized glove. If the measurements do not match, stop here. This glove will not fit.

#### Trimming and Fitting Instructions

1. Apply appropriate sized glove.
2. If the glove fingers are too long, with a pen mark the point along the glove fingers at which they should be trimmed.
3. Remove glove.
4. Trim the glove fingers at all points marked in Step 2. Then, reapply the glove.
5. Do not trim while patient has glove donned.

#### Microfine Toe Cap Measuring Instructions

Measure base of toes and match measurements to appropriate sized Toe Cap. If the measurement does not match, stop here. This Toe Cap will not fit.

#### Trimming and Fitting Instructions

1. Apply appropriate sized toe cap.
2. If the toe cap toes are too long, with a pen mark the point along the toe cap toes at which they should be trimmed.
3. Remove toe cap.
4. Trim the toe cap toes at all points marked in Step 2. Then, reapply the toe cap.
5. Do not trim while patient has garments donned.

### Gauntlet


#### FarrowWrap® Hand Gauntlet Measuring Instructions

Measure hand at base of fingers, fold of thumb, and wrist. Match measurements to the appropriate sized gauntlet. If the measurements do not match, stop here. This gauntlet will not fit.

If the gauntlet thumbhole is too small, it can be trimmed larger **by a fitter**, but **exercise caution not to cut off too much fabric or cut across seams**. If the thumb swells, a glove will be needed to control the swelling.

#### Fitting Instructions

1. Place the gauntlet on a table in front of you. The thumbhole should be closer to the center of the body for the hand on which it will be placed, with the thinner band further away from you.
2. To attach foam, flip the gauntlet over and attach either the 4 or 8 mm foam using the Velcro® tabs on the foam about 1 cm from the thumbhole.
3. Attach the Velcro® pieces to the appropriate bands on the side opposite the foam leaving about half sticking off the edge.




Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowWrap® 4000 Knee High Order Form

Note: All measurements in cm



Please place a numeric value in the white box to indicate items needed.	FarrowWrap® 4000 Legpiece Size Chart							
	XSmall	Qty	Small	Qty	Medium	Qty	Large	Qty
(C) Widest Calf	25 – 38 cm		32 – 48 cm		37 – 56 cm		43 – 67 cm	
(B) Ankle	17 – 23 cm		19 – 28 cm		25 – 37 cm		29 – 43 cm	
(B-D) Regular	30 – 33 cm		30 – 33 cm		30 – 33 cm		30 – 33 cm	
(B-D) Tall	34 – 37 cm		34 – 37 cm		34 – 37 cm		34 – 37 cm	

Color Options	
Tan	<input type="checkbox"/>
Black	<input type="checkbox"/>

*Default color option is Tan. A matching FarrowHybrid® will be included with the FarrowWrap®4000*

### FarrowWrap® 4000 Legpiece

Obtain the following 3 measurements:


- Widest calf circumference (at **C**).
- Least ankle circumference (at **B**) just above malleoli.
- Length of leg following posterior contour from the least ankle to 2 finger widths below knee crease, **B-D**.

### Note

Each FarrowWrap® 4000 legpiece comes with one (1) color matching FarrowHybrid® Foot Compression Liner. **FarrowHybrid® AD Foot Compression** is part compression stocking (over the foot and ankle) and part sock (ankle to knee) and may be used in place of our footpieces with our legpiece. Please be sure to indicate the color of the FarrowWrap® 4000 desired. If no color is selected, the default color will be tan.

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

## FarrowWrap® Trim-To-Fit Knee High Order Form



## FarrowWrap® TTF Knee High

FarrowWrap® TTF Legpiece		Medium	XLarge	FarrowWrap® TTF Footpiece		Medium
(L4) & (L5) Circumferences		20 – 60 cm	45 – 90 cm	Forefoot (X) circumference		20 – 30 cm
(L3) Circumference		20 – 55 cm	45 – 90 cm	(X) Heel to forefoot length		14 – 21 cm
(L2) Circumference		20 – 45 cm	35 – 80 cm	<b>LITE TTF Footpiece</b>	FWLT-T-FT3	Qty
(B) Ankle Circumference		20 – 40 cm	35 – 70 cm	<b>STRONG TTF Footpiece</b>	FWST-T-FT3	Qty
(B-D) Posterior Leg Length		27 – 35 cm	27 – 35 cm	Please note that you may substitute a FarrowHybrid® Foot Compression for a footpiece. See FarrowHybrid® order form for additional details.		
<b>LITE TTF Legpiece</b>	FWLT-T-LT3	Qty				
<b>STRONG TTF Legpiece</b>	FWST-T-LT	Qty	Qty			

FarrowWrap® TTF garments fit a very wide range of sizes and can be trimmed on site for a perfect fit, and retrimmed if the patient’s edema further reduces.

When measuring the circumferences from B to L5, the distance between each circumference will vary for the different fabrics. For **LITE** legpieces, measure the circumference at **point B** and then measure the circumference every **6.5 cm** until measurement **point D** is reached. For **STRONG** legpieces, measure the circumference at **point B** and then measure the circumference every **5 cm** until measurement **point D** is reached.

**FarrowWrap® TTF (Trim-To-Fit) garments require trimming by a professional before being applied by the patient. When a TTF footpiece and legpiece are purchased together, a free pair of liners is included. Silver liners are supplied by default, but you may request Terry Cloth instead.**

# Directions for FarrowWrap® TTF Legpiece Measuring and Fitting


Live web-based in-services are also available. We highly recommend that you schedule an appointment for an in-service.

## Measuring Instructions for LITE and STRONG Legpieces

1. Following the back of the leg, measure from **B** (least ankle) up to **D** (2 finger widths below the back knee crease). The **B-D** length must be between 27–35 cm. If the **B-D** length falls outside this range, the garment will not fit. Contact Farrow Innovations for alternatives.
2. If measuring for a LITE legpiece go to step 3. If measuring for a STRONG legpiece, skip to step 4.
3. For **LITE** TTF legpieces, measure and record the circumference at **B**, then up the leg every 6.5 cm until measurement **point D** is reached.
4. For **STRONG** TTF legpieces, measure and record the circumference at **B**, then up the leg every 5 cm until measurement **point D** is reached.  
\*The circumferences must within the ranges specified for each size garment. If any of the circumferences fall outside these ranges, the garment will not fit.

## Trimming Instructions for LITE and STRONG Legpieces

1. Trim the bottom most band of the TTF legpiece on each side in accordance with the **B** circumference. Distances between lines are representative of 5 cm. If the circumference is between lines, estimate the distance between the 2 most appropriate lines and trim accordingly. See illustration **trimming at 33 cm**.
2. Trim the bands consecutively up the legpiece corresponding to their respective L2 – L5 circumferences.
3. The separate band can be attached, using the Velcro® tab, to the top of the spine (in back) and moved up or down the band below as needed to increase or decrease the height of the legpiece.
4. Attach the Velcro® pieces to the patient's preferred side with **exactly** half extending over the edge of the band. **FIG. 1**
5. Roll the loose end of the Velcro® back onto the same band. **FIG. 2**
6. Instruct the patient not to lose the Velcro®. Replacement Velcro® sets are available for sale through your supplier.



## Fitting Instructions and Important Information

1. LITE and STRONG garments should be applied distal to proximal with ¼ overlap over the opposing bands (**not including Velcro®**). See **FIG. 3**.
2. If the patient has mild to moderate arterial disease and/or peripheral neuropathy, make sure that the prescribed compression garment is appropriate for the patient. If so, for LITE and STRONG garments after the fitter has appropriately applied the bands for the first time, place a mark with a permanent fabric marker a little beyond the applied Velcro® edge so that the patient will have a landmark when applying the garment the next time. Make sure that the patient understands they need to apply the Velcro® with **exactly** half extending over the edge of the band. **FIG. 4**



As long as the patient's leg is an inverted cone shape,  the ankle being smaller than the calf, and the garment is appropriately trimmed and applied, the LITE garments should apply 20-30mmHg and STRONG garments 30-40mmHg of gradient **resting** compression - with **higher working** compression levels provided by short stretch.

Please note that the resting compression will likely increase if the leg tries to swell as a result of the patient's leg being upright (standing or sitting positions). Advise the patient that walking and/or flexing their calf muscle will augment the calf muscle pump and help reduce any edema. If the upper leg is **significantly** larger than the ankle, all legpieces may slip down. If so, consider Farrow GarmentGrip™ as a solution to help prevent slippage. Patients with wider ankles than calves may not receive gradient resting compression levels with this garment. Bandaging by a wound care or lymphedema professional skilled in the art is advised. If the edema reduces significantly the patient will need to be remeasured and the garment(s) retrimmed corresponding to the new circumferences following the Trimming and Fitting Instructions.

If the patient's edema worsens, their health care provider should reevaluate the patient. If continuing with the garment and the patient does not have any arterial disease, peripheral neuropathy, or physical or mental impairments that would prevent them from adjusting or removing the garment, then you could consider applying the bands near maximal stretch. These garments will allow for about a 20% increase in edema; however, keep in mind that the stated compression ranges for LITE and STRONG will be less accurate when applied to a limb larger than that for which the garment was trimmed.

- It is contraindicated to use TTF legpieces on leg circumferences less than 20 cm.
- Do not mix with non Farrow Brand garments as this may increase risk to the patient and is therefore contraindicated.

# Directions for FarrowWrap® TTF Footpiece Measuring and Fitting


Live web-based in-services are also available. We highly recommend that you schedule an appointment for an in-service.

## Measurement Instructions for LITE and STRONG Footpieces

1. Start by measuring the Straight Distance X (length), from the back of the heel straight along the floor to just before the bunions (the forefoot). Straight Distance X must be between 14-21 cm in length. If the patient's Straight Distance X falls outside of this range, stop here. This garment will not fit.
2. Next, measure the circumference at X, just before the bunions. The circumference at X is only guaranteed to fit between 20-30 cm, though it may be possible to fit circumferences up to 39 cm.
3. Lastly, measure the mid-foot circumference. If it is 5 or more centimeters greater than the circumference at X, stop here. This garment will not fit.
  - It is contraindicated to use this garment on foot circumferences less than 19 – 20 cm or on patients with severe arterial disease.
  - Do not mix with non Farrow Brand garments as this may increase risk to the patient and is as such contraindicated.

## Trimming and Fitting Instructions

1. Trim the length of the footpiece to the Straight Distance X on the garment. The first line along the length represents 14 cm, and each consecutive line thereafter represents 1 additional cm.
2. Using the two numbers on either side of the footpiece closest to the toes, count out from that number and trim the width closest to the circumference at X on both sides of the footpiece on or parallel to the lines on each side of the footpiece. The distance between lines represents 2 cm. If the first line were 19, the next would be 21 and so forth. If only 1 cm increment is needed, cut in between and parallel to the lines. **FIG. 1** - trimming for 22 cm circumference.
3. Cut the long strip of Velcro® the same length as the straight edge you cut in step 2. Now, fold it in half and cut it in half. **FIG. 2**
4. To create the flaps for the Velcro®, fold the straight edge of the footpiece, which you just cut, on either side in half with lines facing up. Cut perpendicularly to the edge along the folded crease of the garment 2 cm deep or 3 lines. Do this on both sides. **FIG. 3**
5. Attach the Velcro® to the flaps created on patient's preferred side with half extending past the flap edge. Then fold the Velcro® back over itself onto the garment in preparation for putting it on. Make sure the patient understands that it is important they do the same. **FIG. 4**
6. **For LITE TTF footpieces**, trim ankle bands long enough so that the Velcro® can attach to the opposite side of the foot without touching any underlying Velcro®. Attach the included Velcro® pieces to both sides of the ankle band with half extending past the edge of the band. Then apply. **FIG. 5**
7. **For STRONG TTF footpieces**, to apply additional compression to the foot leave the ankle bands long. The first ankle band should be wrapped across the top of the foot to the opposite side closest to the toes, and can then be wrapped around that portion of the foot again. Affix with the Velcro®. The second ankle band should also be wrapped across to the opposite side of the foot, but closer to the ankle for additional compression to that part of the foot, or it could be used to wrap around the ankle at the discretion of the therapist. See **FIG. 6**. If no additional compression is desirable, they can be trimmed and applied as in step 6.
8. Instruct the patient not to lose the Velcro®. Replacement Velcro® sets are available for sale through your supplier.



**FIG. 1**



**FIG. 2**



**FIG. 3**



**FIG. 4**



**FIG. 5**



**FIG. 6**



**FIG. 7**

## Notes for TTF Footpieces


- When properly trimmed and applied, the LITE garment should provide 20-30 mmHg resting compression, and the STRONG 30-40 mmHg resting compression. (Working compression levels may be higher.)
- If the patient has mild to moderate arterial disease and/or peripheral neuropathy, STRONG is NOT recommended. LITE may be considered, but make sure it is appropriately trimmed and applied, with no more than 1 cm of overlap at the base of the slits. See figure to the right. Once you are sure that the footpiece has been appropriately applied (not too tight), then place a mark with a permanent fabric marker a little beyond the applied Velcro® so that the patient will have a landmark when applying the garment the next time. Make sure the patient understands it is important that the Velcro® is halfway off the edge of the flap. **FIG. 7**
- If the patient's edema reduces significantly, the patient will need to be remeasured and the footpiece retrimmed corresponding to the new circumferences following the steps outlined above. If a mark was placed on the fabric previously as a guide, a new mark placed with a permanent fabric marker will need to be provided. This is especially critical for patients with any arterial disease and/or peripheral neuropathy.
- If the patient's edema worsens, their health care provider should reevaluate the patient. If continuing with the garment, apply the bands near maximal stretch. The garment will allow for some increase in edema; however, keep in mind that the extra stretch will likely increase the compression level, and so great care should be exercised when doing so on patients with arterial disease, peripheral neuropathy, or physical or mental impairment that would prevent them from adjusting the footpiece or communicating discomfort to the care giver.

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowWrap® Trim-To-Fit Arm Order Form



### FarrowWrap® LITE TTF Arm

FarrowWrap® TTF Arm			Medium	XLarge
(G) axilla circumference			25 – 44 cm	30 – 48 cm
(F) widest bicep circumference			25 – 42 cm	25 – 45 cm
(E) elbow crease circumference			20 – 38 cm	25 – 42 cm
(D) widest forearm circumference			20 – 36 cm	20 – 40 cm
(C) wrist crease circumference			15 – 31 cm	15 – 35 cm
(E-G) elbow crease to axilla length			15 – 23 cm	20 – 28 cm
(C-E) wrist to elbow crease circumference			20 – 25 cm	20 – 28 cm
<b>Left TTF Arm</b>	FWLT-T-ATL	Qty		
<b>Right TTF Arm</b>	FWLT-T-ATR	Qty		

**This garment only requires 1 hand to don and doff.**

**FarrowWrap® TTF (Trim-To-Fit) garments require trimming by a professional before being applied by the patient. A free pair of liners is included with each TTF armpiece. Only silver liners are available for this garment.**


If you are a healthcare provider, live web-based in-services are available. We highly recommend that you schedule an appointment for an in-service.

### Measuring Instructions

1. Obtain the following measurements: C-E, E-G (lengths), and circumferences at C, D, E, F, and G and record these on the form. Match measurements to the table on the measurement form. If any of the measurements fall outside of the ranges specified, stop here. The garment will not fit.

### Trimming and Fitting Instructions

1. Trim the length of the garment at the C-E and E-G lines closest to the respective measurement. Once trimmed, the C-E line closest to the wrist will become the C circumference and the E-G line closest to the shoulder will become the G circumference.
2. Next, mark with a pen or pencil a dot on the garment the circumference measurements at C, D, E, F, and G.
3. Cut a **straight line** between each pair of dots you made in step 2. There should be 4 cuts (C-D, D-E, E-F, and F-G). **DO NOT STAIR STEP** the cuts.
4. Next to create the flaps for the Velcro®, from the edge of the garment cut along each line 3.5 cm. The flap closest to wrist could be as narrow as 2.5 cm. The flap closest to the armpit should not be cut any narrower than 5 cm. It may be necessary to combine the second to last flap with the last flap to create a single 7.5 cm wide flap.
5. Apply the correct sized Velcro® pieces for each respective flap leaving half sticking off the edge of the flap.
6. Roll each piece of Velcro® back on itself.
7. Along the opposite side, apply 3 symmetrically shaped pieces of Velcro®, one on the side closest to the arm pit, one on the side closest to the wrist, and one near the center where it will be just a little above the elbow crease. Leave half sticking off the edge. See figure of trimmed garment prepared for donning.
  - It is contraindicated to use this garment on wrist circumferences less than 15 cm or on patients with severe arterial disease.
  - Instruct the patient not to lose the Velcro®. We do sell replacement Velcro® sets.



When properly trimmed and applied and as long as the patient's arm is wider near the shoulder and narrower near the wrist, the patient should receive gradient resting compression between 20-30mmHg. (Working compression levels may be higher.) If the upper arm is significantly larger than the elbow, the garment may slip down. If so, consider Farrow GarmentGrip™ as a solution to help prevent slippage. Patients with wider wrists or forearms than upper arms may not receive gradient resting compression levels with this garment. Bandaging by a wound care or lymphedema professional skilled in the art is advised.

If the patient has wrist circumference between 15-20 cm and has arterial disease that affects the arms, make sure that the prescribed compression garment is appropriate for the patient. If so, after applying the bands with exactly 2.5 cm (1 inch) of overlap at the slit, consider placing a mark with a permanent fabric marker a little beyond the applied Velcro® so that the patient will have a landmark when applying the garment the next time.



Please note that the resting compression may increase if the patient's arm swells while the garment is applied. Advise the patient that moving their arm will help augment the arm's muscle pump and help reduce any edema.

If the patient's edema reduces significantly, the patient will need to be remeasured and the garment retrimmed corresponding to the new circumferences following steps 2-7. If a mark was placed on the fabric previously as a guide, a new mark placed with a permanent fabric marker will need to be provided. This is especially critical for patients with any arterial disease and/or peripheral neuropathy affecting the arms.


If the patient's edema worsens, their health care provider may need to reevaluate the patient. If continuing with the garment, apply the garment near maximal stretch. So applied, the garment will fit a slightly larger arm, but keep in mind that the extra stretch will likely increase the compression level, and so great care should be exercised when doing so on patients with arterial disease and/or peripheral neuropathy affecting the arm.

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowWrap® Liners Order Form



All silver and Terry Cloth liners are sold as pairs.

**Farrow Silver Liners** are issued standard with most garments. Patients with legs > 70 cm or arms > 58 cm in circumference or that have adverse reactions to silver will need to order a Terry Cloth liner instead, which can be substituted free.

Farrow Silver Liners		Sizes			
		Small (Flesh Color)		Large (Silver Color)	
Calf (C) or Arm Circumference (F¹)		15 – 60 cm		52 – 70 cm	
Silver AD Liner	FLS_-AD_	Qty		Qty	
Silver Arm Liner	FLSF-A2	Qty			
Silver Arm Hand Liner	FLSF-AH2	Qty			

**Terry Cloth** is used on patients who are allergic or adversely react to silver, do not fit into a silver liner, who need a little more cushioning than our silver liners, or who desire an open toe. The toe box can also be sewn shut if necessary. There is no extra charge for this but **MUST** be communicated to us in the space below as the default option is open toe.


Terry Cloth Liners		Sizes					
		Small		Medium		Large	
Extremity Circumference		< 40 cm		40 – 70 cm		70 – 120 cm	
Arm (F¹)	FLTG-ARM	Qty		Qty		Qty	
AD Knee High (C)	FLTG-AD	Qty		Qty		Qty	
Are the AD liners a SUBSTITUTE to our standard liners? (if you circle "YES" you will not be charged.)						YES / NO	
AG Thigh High (G)	FLTG-AG	Qty		Qty		Qty	
Closed toe option? (available for knee high and thigh high) (default is open toe. No extra charge for closed toe)						YES / NO	

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowFoam™ Liners Order Form



All FarrowFoam™ liners are sold as eaches.

**FarrowFoam™ Liners** are used under a Farrow brand compression garment. When ordering a FarrowFoam™ Liner you **MUST** add 10 cm to all leg circumferences and 5 cm to foot circumferences when choosing the correct-sized Farrow leg or foot garment to go over it. On the order form, please indicate the actual measurements and then the measurements needed to fit over the FarrowFoam™ liner.

FarrowFoam™ Smooth Foam AD Liner		Sizes					
		XSmall	Small	Medium	Large	XLarge*	XXLarge*
(D) Below Knee		37 – 43 cm	42 – 48 cm	45 – 51 cm	50 – 56 cm	55 – 61 cm	60 – 66 cm
(C) Mid calf		32 – 38 cm	37 – 43 cm	42 – 48 cm	47 – 53 cm	52 – 58 cm	57 – 63 cm
(B) Least Ankle		18 – 23 cm	21 – 26 cm	23 – 28 cm	26 – 31 cm	28 – 33 cm	31 – 36
(A-D) Length		38 cm	39 cm	41 cm	42 cm	43 cm	45 cm
<b>Side Closure AD Liner</b>	FLFM-AD	Qty	Qty	Qty	Qty	Qty	Qty

\*XL and XXL sizes are subject to price increases.



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## Accessories Order Form

<b>GarmentGrip™</b>			
10cm x 60cm Tan (Arm)	FGG-T-PUT010X060	Qty	
14cm x 80cm Tan (Thigh)	FGG-T-NPT014X080	Qty	
14cm x 120cm Tan (Thigh)	FGG-T-NPT014X120	Qty	

<b>Farrow Short-Stretch™</b>			
8cm x 150cm Tan Webbing	FSSWB-08X150	Qty	
10cm x 150cm Strong Band	FSSSTT-10X150	Qty	

<b>Velcro® Packs</b>			
For LITE TTF Armpiece	FWLT-T-AVP	Qty	
For TTF Footpiece	FW-T-FVP	Qty	
For STRONG TTF Legpiece	FWST-T-LVP	Qty	
For LITE TTF Legpiece	FWLT-T-LVP	Qty	
For STRONG TTF Thighpiece	FWST-T-TVP	Qty	
For BASIC OTS Legpiece	FWBA-O-LVP	Qty	
For BASIC OTS Footpiece	FWBA-O-FVP	Qty	

<b>Velcro® &amp; Foam Packs</b>		<b>Size</b>			
		<b>X-Small – Small</b>		<b>Medium – Large</b>	
For LITE Hand Gauntlet	FWLT-O-HGRP	Qty		Qty	