



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

# FarrowWrap® Prescription Form

## Prescribing Information



	Left	Right
Thigh	<input type="checkbox"/>	<input type="checkbox"/>
Leg	<input type="checkbox"/>	<input type="checkbox"/>
Foot	<input type="checkbox"/>	<input type="checkbox"/>
Toes	<input type="checkbox"/>	<input type="checkbox"/>



Mild to Moderate Edema

Moderate to Severe Edema\*

\*Be cautious when prescribing "Moderate to Severe Edema" for patients with PAD.

If mixed severity, please specify in special instructions below.

	Left	Right
Arm	<input type="checkbox"/>	<input type="checkbox"/>
Hand	<input type="checkbox"/>	<input type="checkbox"/>
Fingers	<input type="checkbox"/>	<input type="checkbox"/>



Special Instructions: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_


Related Diagnoses: \_\_\_\_\_

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowHybrid® Foot Compression Order Form



**FarrowHybrid® AD Foot Compression** is part compression stocking (over the foot and ankle) and part sock (ankle to knee) and may be used in place of our footpieces with our legpiece. Easy to apply and remove, they are comfortable and readily fit in patients' shoes.

FarrowHybrid® AD		SKU	Small	Medium	Large	XLarge
Male Shoe Size			3 – 6	7 – 9	10 – 12	13 – 15
Female Shoe Size			4 – 7	8 – 10	11 – 13	14 – 16
<b>FarrowHybrid® ADI</b>	FLFH-ADI	Qty	Qty	Qty	Qty	
<b>FarrowHybrid® ADI Wide</b>	FLFH-ADIW	Qty	Qty	Qty	Qty	
<b>FarrowHybrid® ADIIS</b>	FLFH-ADIIS	Qty	Qty	Qty	Qty	
<b>FarrowHybrid® ADIIS Wide</b>	FLFH-ADIISW	Qty	Qty	Qty	Qty	

Color Options	
Taupe	
Black	

*Black available only in FarrowHybrid ADI Wide. If no color is selected, Taupe is the default option.*

### FarrowHybrid® Foot Compression Size Selector

	20 – 30mmHg	30 – 40mmHg
A¹ Circumference 20 – 27cm (standard)	FLFH-ADI	FLFH-ADIIS
A¹ Circumference 28 – 35cm (wide)	FLFH-ADIW	FLFH-ADIISW

FarrowHybrids cover a wide range of circumferences in standard and wide sizes. The maximum leg circumference must be less than 60cm. FarrowHybrid® ADI and ADIW provide Class I (20-30mmHg) compression and FarrowHybrid® ADIIS and ADIISW provide Class II (30-40mmHg) compression. The FarrowHybrid® ADIIS foot compression CONTAINS SILVER and should not be used on patients allergic to silver.

### FarrowHybrid® AD Foot Compression