



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

# FarrowWrap® Prescription Form

## Patient Information



	Left	Right
Thigh	<input type="checkbox"/>	<input type="checkbox"/>
Leg	<input type="checkbox"/>	<input type="checkbox"/>
Foot	<input type="checkbox"/>	<input type="checkbox"/>
Toes	<input type="checkbox"/>	<input type="checkbox"/>

<input type="checkbox"/>	Mild to Moderate Edema
<input type="checkbox"/>	Moderate to Severe Edema*

\*Be cautious when prescribing "Moderate to Severe Edema" for patients with PAD.

If mixed severity, please specify in special instructions below.

	Left	Right
Arm	<input type="checkbox"/>	<input type="checkbox"/>
Hand	<input type="checkbox"/>	<input type="checkbox"/>
Fingers	<input type="checkbox"/>	<input type="checkbox"/>



Special Instructions: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_


Related Diagnoses: \_\_\_\_\_

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## FarrowWrap® Off-The-Shelf Arm



FarrowWrap® OTS Arm		Small	Medium	Large
(G) axilla circumference		22 – 31 cm	29 – 39 cm	32 – 45 cm
(E) elbow crease circumference		20 – 27 cm	25 – 34 cm	30 – 40 cm
(C) wrist crease circumference		14 – 18 cm	16 – 21 cm	19 – 25 cm
(C-G) Short	<b>40 – 43 cm</b>	Qty		
(C-G) Regular	<b>44 – 47 cm</b>	Qty		
(C-G) Long	<b>48 – 51 cm</b>	Qty		

Options	
<b>SIDE</b>	✓
Left	
Right	
<b>Color</b>	✓
Tan	
Black	

**Left or Right must be selected in the “Options” box.  
Liners will not be issued with this garment.**

*If no color is selected, Tan is the default option.*

### Directions for the Off-The-Shelf Arm

For length, measure from the outside flexor surface of the wrist to the axilla with a slight bend in the arm. Measure circumference at the wrist, elbow crease and axilla. Compare measurements to the size chart to determine correct size. All circumference measurements must correspond to one off-the-shelf size in order for the garment to fit. This garment can be donned and doffed with one hand. **Liners will not be issued with this garment. The OTS Arm is available in LITE fabric only.**

Other: \_\_\_\_\_