



Absolute Medical. 1843 W Hubbard St #2A
Chicago, IL 60622
Tel: 312 233 2207 Fax: 866 860 9358
To Order Online: orders@absolutemedical.com
Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____


Fitter Last Name: _____ Fitter First Name: _____

Date: _____

ADVI Ready Made Order Form

ORDER SPECIFICATIONS:					
<input type="checkbox"/> Quote Only	<input type="checkbox"/> Quote & Proceed	<input type="checkbox"/> Reseller Pricing	<input type="checkbox"/> Retail Pricing		
PREFERRED SHIPPING METHOD:					
<input type="checkbox"/> UPS 3 rd Party Billing (3 rd Party Billing will have a \$5 handling fee)					
<input type="checkbox"/> Most Cost Effective	<input type="checkbox"/> USPS Priority Mail® Small (\$6) or Medium (\$12) Flat Rate boxes	<input type="checkbox"/> UPS® Ground (\$13*)	<input type="checkbox"/> UPS 3 Day Select® (\$26*) <i>Not available for AK or HI</i>	<input type="checkbox"/> UPS 2nd Day Air® (\$30*)	<input type="checkbox"/> UPS Next Day Air® (\$57*)

* UPS shipping prices listed are for contiguous 48 US states only. Shipments to residential addresses will have additional \$7 charge.



Size	
<input type="checkbox"/> XSmall	<input type="checkbox"/> XSmall MAX
<input type="checkbox"/> Small	<input type="checkbox"/> Small MAX
<input type="checkbox"/> Medium	<input type="checkbox"/> Medium MAX
<input type="checkbox"/> Large	<input type="checkbox"/> Large MAX
<input type="checkbox"/> XLarge	<input type="checkbox"/> XLarge MAX
<input type="checkbox"/> 2XLarge	<input type="checkbox"/> 2XLarge MAX

Length
<input type="checkbox"/> XShort
<input type="checkbox"/> Short
<input type="checkbox"/> Regular
<input type="checkbox"/> Long
<input type="checkbox"/> XLong

Side
<input type="checkbox"/> Right
<input type="checkbox"/> Left

Fabric/Color	
Organic Cotton/Lycra®	<input type="checkbox"/> Ivory
Foam	
Low ILD - 2 Blend only	
JoViJacket	
Nylon/Spandex	<input type="checkbox"/> White
Additional Charge Options	
<input type="checkbox"/> 8cm x 5m Lenkelast® Medium Stretch Bandage	
<input type="checkbox"/> Prepaid Reduction Option	

Comments/Quantity:



Absolute Medical. 1843 W Hubbard St #2A
Chicago, IL 60622
Tel: 312 233 2207 Fax: 866 860 9358
To Order Online: orders@absolutemedical.com
Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

ADVI Ready Made Order Form

FOR FITTER/THERAPIST USE ONLY

ADVI SIZE CHART			
Size	Centimeters		
	Least Ankle Circumference	Widest Calf Circumference	Least Knee Circumference
	"B" Landmark	"C" Landmark	"D" Landmark
XSmall	19 - 20.5	33 - 35.5	29 - 32
XSmall MAX	24 - 25.5	35.5 - 38	30.5 - 33
Small	20.5 - 22	35.5 - 38	32 - 34.5
Small MAX	25.5 - 27	38 - 40.5	33 - 35.5
Medium	22 - 23	38 - 40.5	34.5 - 37
Medium MAX	27 - 28	40.5 - 43	35.5 - 38
Large	23 - 24	40.5 - 43	37 - 39.5
Large MAX	28 - 29	43 - 45.5	38 - 41
XLarge	24 - 25.5	43 - 45.5	39.5 - 42
XLarge MAX	29 - 30.5	45.5 - 48	41 - 43
2XLarge	25.5 - 27	45.5 - 48	42 - 44.5
2XLarge MAX	30.5 - 32	48 - 51	43 - 46
MAX sizes are 5 cm (2 in) larger at the ankle, 2.5 cm (1 in) larger at the calf, 1.3 cm (1/2 in) larger at the knee.			

ADVI LEG LENGTH CHART	
Floor to Tibial Tuberosity (A to D)	
XShort	34 cm
Short	38
Regular	42
Long	46
XLong	50