



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____


Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Classic Arm Sleeve Ready Made Order Form

ORDER SPECIFICATIONS:					
<input type="checkbox"/> Quote Only	<input type="checkbox"/> Quote & Proceed	<input type="checkbox"/> Reseller Pricing	<input type="checkbox"/> Retail Pricing		
PREFERRED SHIPPING METHOD:					
<input type="checkbox"/> UPS 3 rd Party Billing (3 rd Party Billing will have a \$5 handling fee)					
<input type="checkbox"/> Most Cost Effective	<input type="checkbox"/> USPS Priority Mail [®] Small (\$6) or Medium (\$12) Flat Rate boxes	<input type="checkbox"/> UPS [®] Ground (\$13*)	<input type="checkbox"/> UPS 3 Day Select [®] (\$26*) <i>Not available for AK or HI</i>	<input type="checkbox"/> UPS 2nd Day Air [®] (\$30*)	<input type="checkbox"/> UPS Next Day Air [®] (\$57*)

* UPS shipping prices listed are for contiguous 48 US states only. Shipments to residential addresses will have additional \$7 charge.



Style	Size		Length	Fabric/Color	
<input type="checkbox"/> Standard	<input type="checkbox"/> XSmall	<input type="checkbox"/> XSmall MAX	<input type="checkbox"/> XShort	Polartec [®] Power Dry [®]	
<input type="checkbox"/> Slimline	<input type="checkbox"/> Small	<input type="checkbox"/> Small MAX	<input type="checkbox"/> Short	<input type="checkbox"/> Royal Blue	
Side			<input type="checkbox"/> Regular	<input type="checkbox"/> Pink	
<input type="checkbox"/> Right			<input type="checkbox"/> Long	<input type="checkbox"/> Black	
<input type="checkbox"/> Left			<input type="checkbox"/> XLong	Organic Cotton/Lycra [®]	
			JoViJacket		
			Nylon/Spandex		
			<input type="checkbox"/> Black		
			<input type="checkbox"/> White		

Additional Charge Options
<input type="checkbox"/> Dycem [®] - donning aid
<input type="checkbox"/> Easy Slide - donning aid
<input type="checkbox"/> Prepaid Reduction Option

Comments/Quantity:



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Classic Arm Sleeve Ready Made Order Form

FOR FITTER/THERAPIST USE ONLY

CLASSIC ARM SIZE CHART			
Size	<i>Centimeters</i>		
	Least Wrist Circumference	Least Elbow Circumference	Axilla Circumference
	<i>"C"</i> Landmark	<i>"E"</i> Landmark	<i>"G"</i> Landmark
XSmall	14 - 15.5	20 - 23	23 - 27.5
XSmall MAX		22 - 25	27 - 31
Small	15.5 - 16.5	23 - 26	27.5 - 32
Small MAX		25 - 28	31 - 36
Medium	16.5 - 18	26 - 29	32 - 36
Medium MAX		28 - 31	36 - 40
Large	18 - 19	29 - 33	36 - 41
Large MAX		31 - 35	40 - 44.5
XLarge	19 - 20.5	33 - 37.5	41 - 45
XLarge MAX		35 - 39.5	44.5 - 49
<i>MAX sizes are 2 cm (3/4 in) larger at the elbow and 4 cm (1 1/2 in) larger at the axilla.</i>			

CLASSIC ARM LENGTH CHART	
Medial Inseam Wrist to Axilla	
XShort	38cm
Short	40.5
Regular	43
Long	45.5
XLong	48