



Absolute Medical, 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Date: _____

COMPREFLEX REDUCE LOWER LEG

PRODUCT INFORMATION

COMPREFLEX REDUCE BK	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
COMPREFLEX REDUCE BOOT	Size		Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

MEASURING INSTRUCTIONS

COMPREFLEX REDUCE BELOW KNEE

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the **Lateral Malleolus** (outer ankle bone), continue to the 25cm dot. Record in the corresponding lines A and C (Left/Right). Measure length from the **Lateral Malleolus** (outer ankle bone) to lateral aspect of the **Fibular Head**, and record in box G.

COMPREFLEX REDUCE BOOT

Measure circumference of ankle bend and heel and record on line I. Measure circumference across metatarsal heads and record on line J. Measure length on medial side of foot from heel to 1st metatarsal head and record on line K.

COMPREFLEX REDUCE BK SIZING CHART

	RANGE 1	RANGE 3
C	<65cm	<90cm
A	<55cm	<75cm
BLACK		
REGULAR	2401-BKR	2403-BKR
TALL	2401-BKT	2403-BKT
BEIGE		
REGULAR	2411-BKR	2413-BKR
TALL	2411-BKT	2413-BKT

COMPREFLEX REDUCE BOOT SIZING CHART

	RANGE 1	RANGE 3
I	30-44cm	39-54cm
J	24-33cm	30-43cm
K	16-25cm	20-27cm
BLACK	2801-BTL	2803-BTL
BEIGE	2811-BTL	2813-BTL

CIRCUMFERENCE

