



Absolute Medical, 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

# COMPREFLEX REDUCE UPPER LEG

## PRODUCT INFORMATION

COMPREFLEX REDUCE TC	Size	Length	Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE
COMPREFLEX REDUCE KNEE	Size		Item #	Color
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				<input type="checkbox"/> BLACK <input type="checkbox"/> BEIGE

## MEASURING INSTRUCTIONS


### COMPREFLEX REDUCE THIGH COMPONENT

Measure circumferences where indicated by the black dots, beginning with the dot 5cm above the lateral aspect **Top of Patella** and continue to the 23cm dot. Record in the corresponding lines D and F (Left/Right). Measure length from the lateral aspect **Top of Patella** to lateral aspect **Gluteal Fold**, and record in box G1.

### COMPREFLEX REDUCE KNEE

Measure circumferences where indicated by the black dots, beginning with the dot 5cm below the lateral aspect **Fibular Head**. Continue with the dot 5cm above the lateral aspect **Top of Patella**. Record in the corresponding lines C1 and D (Left/Right).

## CIRCUMFERENCE



### COMPREFLEX REDUCE TC SIZING CHART

	RANGE 1	RANGE 3
<b>F</b>	<80cm	<125cm
<b>D</b>	<65cm	<110cm

BLACK	LEFT		
	REGULAR	2401-TCR-L	2403-TCR-L
TALL	2401-TCT-L	2403-TCT-L	
BLACK	RIGHT		
	REGULAR	2401-TCR-R	2403-TCR-R
TALL	2401-TCT-R	2403-TCT-R	
BEIGE	LEFT		
	REGULAR	2411-TCR-L	2413-TCR-L
TALL	2411-TCT-L	2413-TCT-L	
BEIGE	RIGHT		
	REGULAR	2411-TCR-R	2413-TCR-R
TALL	2411-TCT-R	2413-TCT-R	

### COMPREFLEX REDUCE KNEE SIZING CHART

	RANGE 1	RANGE 3
<b>D</b>	<75cm	<115cm
<b>C1</b>	<65cm	<90cm
BLACK	2401-KP	2403-KP
BEIGE	2411-KP	2413-KP