

Patient Last Name: \_\_\_\_\_

Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_

Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_



Absolute Medical. 1843 W Hubbard St #2A  
Chicago, IL 60622  
Tel. 312 233 2207 Fax 866 860 9358  
To Order Online: orders@absolutemedical.net  
Our website: absolutemedical.com

# Lower Extremity Order Form

Quantity/Class	CCL 1 18-21 mmHg*	CCL 2 23-32 mmHg*	CCL 3 34-46 mmHg*	CCL 3F 34-46 mmHg*	CCL 4 49-70 mmHg*	CCL 4S 60-90 mmHg*
Left						
Right						
Pressure panty Elvarex® (Body Bandage)						

Form 57021 must accompany this form.		
Circum. (c)	Length (l)	Length (l)
cT	K2-T	lT
cH	K1-T	lH

### Basic styles

- AD Knee
- AF Mid-Thigh
- AG Thigh
- AG-T (AG with chapstyle)
  - Piece
  - Pair
- AG-HT 1½ Leg Pantyhose
- AT Pantyhose

### Color

- Beige
- Black
- Cocoa
- Navy
- Grey
- Cranberry
- Graphite†
- Henna†
- Denim†

### Special options

- SoftFit band** (Only Available CCL 1-3, AD only)
- Silicone dotted band 2.5cm:**
  - Top
  - Inside
  - Inside ¾
- Silicone dotted band 5cm:**
  - Top
  - 5cm Top waistband
  - Inside
  - Inside ¾
- Ribbed Elastic Band 5cm:**
  - On Top, Waist High Only
- Sensitive Band 5cm:**
  - On Top, Waist High Only
- Vertical Silicone Strips AG\*\*:**
  - Front
  - Back
  - Both
- Zipper B to D only:**
  - Inside (medial)
  - Outside (Lateral)
- Zipper E to G only:**
  - Inside (medial)
  - Outside (Lateral)
- Knee Comfort Zone** (Not available in CCL 1)
- Pocket**
  - Instep (Not Available with Profile)
  - Back of Knee\*\*
  - All four sides closed
- Fly for Men**
- Open Pubis**
- Adjustable Waistband**
- Top Comfort Zone**

### Variations

- B'G
- BG
- FG (leg extension)
- GT Biker Shorts
- B1G-T Capri


### Options

- T-Heel (Class 2-3 forte only)
- Profile (Not available with pocket instep)

### Profile Measuring Points

Point P1: horizontal circumference from the inner to the outer malleolus.

Point P2: "stirrup" or vertical circumference from the inner to the outer malleolus.





† Cannot combine with a different seam color.

### Comments:

Circumference (c)		Length (l): Taken from each landmark to floor	
Left	Right	Left	Right
cG		lG	
cF		lF	
cE		lE	
cD		lD	
cC		lC	
cB1		lB1	
cB		lB	
cY		lA (medial)	
cA		lA (lateral)	

- Straight Open Toe**  
Lateral Length \_\_\_\_\_ cm
- Slant Open Toe**  
Medial Length \_\_\_\_\_ cm
- Slant Closed Toe**  
Medial Length \_\_\_\_\_ cm
- Straight Closed Toe**  
Total Foot Length \_\_\_\_\_ cm
- Lateral Length \_\_\_\_\_ cm
- Lateral Length \_\_\_\_\_ cm
- Total Foot Length \_\_\_\_\_ cm



**CAUTION:** This product contains natural rubber latex which may cause allergic reactions.

\* Design Pressure  
\*\* Not available with Knee Comfort Zone