



ComfyBoot

Ready Made Order Form

FOR FITTER/THERAPIST USE ONLY

COMFYBOOT SIZE CHART			
Size	Centimeters		
	Least Ankle Circumference	Widest Calf Circumference	Least Knee Circumference
	"B" Landmark	"C" Landmark	"D" Landmark
XSmall	19 - 20.5	33 - 35.5	29 - 32
XSmall MAX	24 - 25.5	35.5 - 38	30.5 - 33
Small	20.5 - 22	35.5 - 38	32 - 34.5
Small MAX	25.5 - 27	38 - 40.5	33 - 35.5
Medium	22 - 23	38 - 40.5	34.5 - 37
Medium MAX	27 - 28	40.5 - 43	35.5 - 38
Large	23 - 24	40.5 - 43	37 - 39.5
Large MAX	28 - 29	43 - 45.5	38 - 41
XLarge	24 - 25.5	43 - 45.5	39.5 - 42
XLarge MAX	29 - 30.5	45.5 - 48	41 - 43
2XLarge	25.5 - 27	45.5 - 48	42 - 44.5
2XLarge MAX	30.5 - 32	48 - 51	43 - 46
MAX sizes are 5 cm (2 in) larger at the ankle, 2.5 cm (1 in) larger at the calf, 1.3 cm (1/2 in) larger at the knee.			

LEG LENGTH CHART	
Floor to Tibial Tuberosity (A to D)	
XShort	34 cm
Short	38
Regular	42
Long	46
XLong	50

FOOT LENGTH CHART	
Length is back of heel to tip of longest toe	
Short	21.6-24.3cm
Regular	24.3-26.6
Long	26.6-29.2