



Absolute Medical. 1843 W Hubbard St #2A  
Chicago, IL 60622  
Tel: 312 233 2207 Fax: 866 860 9358  
To Order Online: [orders@absolutemedical.com](mailto:orders@absolutemedical.com)  
Our website: [absolutemedical.com](http://absolutemedical.com)


Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
Date: \_\_\_\_\_

## ComfyBoot

### Ready Made Order Form

ORDER SPECIFICATIONS:					
<input type="checkbox"/> Quote Only	<input type="checkbox"/> Quote & Proceed	<input type="checkbox"/> Reseller Pricing	<input type="checkbox"/> Retail Pricing		
PREFERRED SHIPPING METHOD:					
<input type="checkbox"/> UPS 3 <sup>rd</sup> Party Billing (3 <sup>rd</sup> Party Billing will have a \$5 handling fee)					
<input type="checkbox"/> Most Cost Effective	<input type="checkbox"/> USPS Priority Mail® Small (\$6) or Medium (\$12) Flat Rate boxes	<input type="checkbox"/> UPS® Ground (\$13*)	<input type="checkbox"/> UPS 3 Day Select® (\$26*) <i>Not available for AK or HI</i>	<input type="checkbox"/> UPS 2nd Day Air® (\$30*)	<input type="checkbox"/> UPS Next Day Air® (\$57*)

\* UPS shipping prices listed are for contiguous 48 US states only. Shipments to residential addresses will have additional \$7 charge.



Size	
<input type="checkbox"/> XSmall	<input type="checkbox"/> XSmall MAX
<input type="checkbox"/> Small	<input type="checkbox"/> Small MAX
<input type="checkbox"/> Medium	<input type="checkbox"/> Medium MAX
<input type="checkbox"/> Large	<input type="checkbox"/> Large MAX
<input type="checkbox"/> XLarge	<input type="checkbox"/> XLarge MAX
<input type="checkbox"/> 2XLarge	<input type="checkbox"/> 2XLarge MAX

Length
<input type="checkbox"/> XShort
<input type="checkbox"/> Short
<input type="checkbox"/> Regular
<input type="checkbox"/> Long
<input type="checkbox"/> XLong

Foot Length
<input type="checkbox"/> Short
<input type="checkbox"/> Regular
<input type="checkbox"/> Long

Fabric/Color	
Organic Cotton/ Lycra® and Nylon/ Spandex Techsheen	<input type="checkbox"/> Black/Black <input type="checkbox"/> Ivory/Buff

Additional Charge Options
<input type="checkbox"/> Dycem® - donning aid
<input type="checkbox"/> Prepaid Reduction Option

Comments/Quantity:

---



---



---



---



---



---



# ComfyBoot

## Ready Made Order Form

FOR FITTER/THERAPIST USE ONLY

COMFYBOOT SIZE CHART			
Size	Centimeters		
	Least Ankle Circumference	Widest Calf Circumference	Least Knee Circumference
	"B" Landmark	"C" Landmark	"D" Landmark
XSmall	19 - 20.5	33 - 35.5	29 - 32
XSmall MAX	24 - 25.5	35.5 - 38	30.5 - 33
Small	20.5 - 22	35.5 - 38	32 - 34.5
Small MAX	25.5 - 27	38 - 40.5	33 - 35.5
Medium	22 - 23	38 - 40.5	34.5 - 37
Medium MAX	27 - 28	40.5 - 43	35.5 - 38
Large	23 - 24	40.5 - 43	37 - 39.5
Large MAX	28 - 29	43 - 45.5	38 - 41
XLarge	24 - 25.5	43 - 45.5	39.5 - 42
XLarge MAX	29 - 30.5	45.5 - 48	41 - 43
2XLarge	25.5 - 27	45.5 - 48	42 - 44.5
2XLarge MAX	30.5 - 32	48 - 51	43 - 46
MAX sizes are 5 cm (2 in) larger at the ankle, 2.5 cm (1 in) larger at the calf, 1.3 cm (1/2 in) larger at the knee.			

LEG LENGTH CHART	
Floor to Tibial Tuberosity (A to D)	
XShort	34 cm
Short	38
Regular	42
Long	46
XLong	50

FOOT LENGTH CHART	
Length is back of heel to tip of longest toe	
Short	21.6-24.3cm
Regular	24.3-26.6
Long	26.6-29.2