





## InnaBoot - Lower Leg Ready Made Order Form

FOR FITTER/THERAPIST USE ONLY

<b>InnaBoot Size Chart</b>			
Size	Least Ankle Circumference	Widest Calf Circumference	Least Knee Circumference
	“B” Landmark	“C” Landmark	“D” Landmark
<b>XSmall</b>	19 - 20.5cm	33 - 35.5cm	29 - 35.5cm
<b>XSmall MAX</b>	24 - 25.5	35.5 - 38	30.5 - 36.8
<b>Small</b>	20.5 - 22	35.5 - 38	32 - 38
<b>Small MAX</b>	25.5 - 27	38 - 40.5	33 - 39.3
<b>Medium</b>	22 - 23	38 - 40.5	34.5 - 40.6
<b>Medium MAX</b>	27 - 28	40.5 - 43	35.5 - 41.9
<b>Large</b>	23 - 24	40.5 - 43	37 - 43.1
<b>Large MAX</b>	28 - 29	43 - 45.5	38 - 44.4
<b>XLarge</b>	24 - 25.5	43 - 45.5	39.5 - 45.7
<b>XLarge MAX</b>	29 - 30.5	45.5 - 48	41 - 46.9
<b>2XLarge</b>	25.5 - 27	45.5 - 48	42 - 48.2
<b>2XLarge MAX</b>	30.5 - 32	48 - 51	43 - 49.5
MAX Sizes are 5 cm larger at the ankle, 2.5 cm larger at the calf, 1.3 cm larger at the knee.			

<b>InnaBoot Leg Length Chart</b>	
Medial Inseam (Floor to Least Knee)	
<b>Short</b>	38 cm
<b>Regular</b>	42
<b>Long</b>	46
<b>XLong</b>	50

<b>InnaBoot Foot Length Chart</b>	
Length is back of heel to tip of longest toe	
<b>Short</b>	21.6-24.3cm
<b>Regular</b>	24.3-26.6
<b>Long</b>	26.6-29.2