



Absolute Medical. 1843 W Hubbard St #2A
Chicago, IL 60622
Tel: 312 233 2207 Fax: 866 860 9358
To Order Online: orders@absolutemedical.com
Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____
Fitter Last Name: _____ Fitter First Name: _____
Date: _____

Transitions - Bandage/Pump Liners

Ready Made Order Form

ORDER SPECIFICATIONS:					
<input type="checkbox"/> Quote Only	<input type="checkbox"/> Quote & Proceed	<input type="checkbox"/> Reseller Pricing	<input type="checkbox"/> Retail Pricing		
PREFERRED SHIPPING METHOD:					
<input type="checkbox"/> UPS 3 rd Party Billing (3 rd Party Billing will have a \$5 handling fee)					
<input type="checkbox"/> Most Cost Effective	<input type="checkbox"/> USPS Priority Mail® Small (\$6) or Medium (\$12) Flat Rate boxes	<input type="checkbox"/> UPS® Ground (\$13*)	<input type="checkbox"/> UPS 3 Day Select® (\$26*) <i>Not available for AK or HI</i>	<input type="checkbox"/> UPS 2nd Day Air® (\$30*)	<input type="checkbox"/> UPS Next Day Air® (\$57*)

* UPS shipping prices listed are for contiguous 48 US states only. Shipments to residential addresses will have additional \$7 charge.

Transitions garments are made of Ivory Organic Cotton/Lycra®

Full Arm	
Size	Length
<input type="checkbox"/> Small	<input type="checkbox"/> Short
<input type="checkbox"/> Medium	<input type="checkbox"/> Regular
<input type="checkbox"/> Large	<input type="checkbox"/> Long
<input type="checkbox"/> XLarge	<input type="checkbox"/> XLong

Lower Leg	
Size	Length
<input type="checkbox"/> Small	<input type="checkbox"/> Short
<input type="checkbox"/> Medium	<input type="checkbox"/> Regular
<input type="checkbox"/> Large	<input type="checkbox"/> Long
<input type="checkbox"/> XLarge	<input type="checkbox"/> XLong
<input type="checkbox"/> 2XLarge	

Upper Leg	
Size	Length
<input type="checkbox"/> Small	<input type="checkbox"/> Short
<input type="checkbox"/> Medium	<input type="checkbox"/> Regular
<input type="checkbox"/> Large	<input type="checkbox"/> Long
<input type="checkbox"/> XLarge	<input type="checkbox"/> XLong
<input type="checkbox"/> 2XLarge	

Full Leg	
Size	Length
<input type="checkbox"/> Small	<input type="checkbox"/> Short
<input type="checkbox"/> Medium	<input type="checkbox"/> Regular
<input type="checkbox"/> Large	<input type="checkbox"/> Long
<input type="checkbox"/> XLarge	<input type="checkbox"/> XLong
<input type="checkbox"/> 2XLarge	

Comments/Quantity:



Absolute Medical. 1843 W Hubbard St #2A
Chicago, IL 60622
Tel: 312 233 2207 Fax: 866 860 9358
To Order Online: orders@absolutemedical.com
Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Transitions - Bandage/Pump Liners

Ready Made Order Form

FOR FITTER/THERAPIST USE ONLY



FULL ARM SIZE CHART	
Size	Widest Bicep Circumference
Small	23 - 30.5 cm
Medium	30.5 - 38
Large	38 - 46
XLarge	46 - 56

WRIST TO AXILLA LENGTH	
Size	Centimeters
Short	40.5
Regular	43
Long	45.5
XLong	48



LOWER LEG SIZE CHART	
Size	Widest Calf Circumference
Small	33 - 38 cm
Medium	38 - 46
Large	46 - 58
XLarge	58 - 66
2XLarge	66 - 74

FLOOR TO KNEE LENGTH	
Size	Centimeters
Short	35
Regular	40
Long	45
XLong	50



UPPER LEG SIZE CHART	
Size	Widest Thigh Circumference
Small	48 - 58 cm
Medium	58 - 74
Large	74 - 89
XLarge	89 - 102
2XLarge	102 - 115

KNEE TO GROIN LENGTH	
Size	Centimeters
Short	32
Regular	36
Long	40
XLong	44



FULL LEG SIZE CHART		
Size	Widest Calf Circumference	Widest Thigh Circumference
Small	33 - 38 cm	48 - 58 cm
Medium	38 - 46	58 - 74
Large	46 - 58	74 - 89
XLarge	58 - 66	89 - 102
2XLarge	66 - 74	102 - 115

FLOOR TO GROIN LENGTH	
Size	Centimeters
Short	69
Regular	74
Long	79
XLong	84