





Absolute Medical. 1843 W Hubbard St #2A
Chicago, IL 60622
Tel: 312 233 2207 Fax: 866 860 9358
To Order Online: orders@absolutemedical.com
Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Jobst Bella Strong Gauntlet Sizing Chart

	 Wrist		 Palm	
Size	Min.	Max.	Min.	Max.
1	5 1/2" (14 cm)	7" (18 cm)	6" (15 cm)	6 3/4" (17 cm)
2	5 1/2" (14 cm)	7" (18 cm)	6 3/4" (17 cm)	7 1/2" (19 cm)
3	7" (18 cm)	8 1/4" (21 cm)	6 3/4" (17 cm)	7 1/2" (19 cm)
4	5 1/2" (14 cm)	7" (18 cm)	7 1/2" (19 cm)	8 1/4" (21 cm)
5	7" (18 cm)	8 1/4" (21 cm)	7 1/2" (19 cm)	8 1/4" (21 cm)
6	6 1/4" (16 cm)	7 1/2" (19 cm)	8 1/4" (21 cm)	9" (23 cm)
7	7 1/2" (19 cm)	8 1/2" (22 cm)	8 1/4" (21 cm)	9" (23 cm)
8	7 1/2" (19 cm)	8 1/2" (22 cm)	9" (23 cm)	10 1/4" (26 cm)