



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com


Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

3 Easy Measuring Steps:

1. Length (B-D) – Ankle to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) – Calf (widest)
3. Circumference (B) – Ankle (smallest part)



JOBST® FarrowWrap® 4000 Size Chart

	XS	S	M	L
C Mid-Calf	25-38 cm	32-48 cm	37-56 cm	43-67 cm
B Ankle	17-23 cm	19-28 cm	25-37 cm	29-43 cm
B-D Regular	30-33 cm	30-33 cm	30-33 cm	30-33 cm
B-D Tall	34-37 cm	34-37 cm	34-37 cm	34-37 cm

JOBST® FarrowWrap® 4000

Ref No.	Size	Length	Color	Quantity
7666100	XS	Regular	Tan	
7666101	S	Regular	Tan	
7666102	M	Regular	Tan	
7666103	L	Regular	Tan	
7666108	XS	Regular	Black	
7666109	S	Regular	Black	
7666110	M	Regular	Black	
7666111	L	Regular	Black	
7666104	XS	Tall	Tan	
7666105	S	Tall	Tan	
7666106	M	Tall	Tan	
7666107	L	Tall	Tan	
7666112	XS	Tall	Black	
7666113	S	Tall	Black	
7666114	M	Tall	Black	
7666115	L	Tall	Black	



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

3 Easy Measuring Steps:

1. Length (A-D) – Floor to 2 finger widths below crease of knee (follow contour of leg)
2. Circumference (C) – Calf (widest)
3. Circumference (B) – Ankle (smallest part)

JOBST® FarrowWrap® Basic (30-40 mmHg*)

Ref No.	Size	Length	Quantity
7666000	XS	Regular	
7666001	S	Regular	
7666002	M	Regular	
7666003	L	Regular	
7666004	XL	Regular	
7666005	XS	Tall	
7666006	S	Tall	
7666007	M	Tall	
7666008	L	Tall	
7666009	XL	Tall	

JOBST® FarrowWrap® Basic Size Chart

	XS	S	M	L	XL
C Mid-Calf	36-43 cm	42-50 cm	48-58 cm	53-63 cm	58-63 cm
B Ankle	21-25 cm	25-30 cm	30-36 cm	36-42 cm	42-50 cm
A-D Regular	33-37 cm	35-39 cm	37-41 cm	39-43 cm	39-43 cm
A-D Tall	38-41 cm	40-43 cm	42-45	44-47 cm	44- cm

* The mean compression for an average size ankle
 Please refer to the product label and / or package insert for full instructions on the safe use of these products.