




Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

## FarrowWrap® Liner/Farrowhybrid® Liner Order Form



### JOBST FARROWWRAP® TG® SOFT KNEE - HIGH LINERS

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7666900	<40 cm (C)	TG SOFT AD (LEG) SIZE SMALL	
7666902	40-70 cm (C)	TG SOFT AD (LEG) SIZE MEDIUM	
7666904	70-120 cm (C)	TG SOFT AD (LEG) SIZE LARGE	

### JOBST FARROWWRAP® TG® SOFT THIGH - HIGH LINERS

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7667000	<40 cm (G)	TG SOFT AG (THIGH) SIZE SMALL	
7667002	40-70 cm (G)	TG SOFT AG (THIGH) SIZE MEDIUM	
7667004	70-120 cm (G)	TG SOFT AG (THIGH) SIZE LARGE	

### JOBST® FARROWFOAM™ LINERS

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7666500	X SMALL	SIDE CLOSURE AD LINERS	
7666501	SMALL	SIDE CLOSURE AD LINERS	
7666502	MEDIUM	SIDE CLOSURE AD LINERS	
7666503	LARGE	SIDE CLOSURE AD LINERS	
7666504	X LARGE	SIDE CLOSURE AD LINERS	
7666505	XX LARGE	SIDE CLOSURE AD LINERS	

### JOBST FARROWFOAM™ LINERS SIZING CHART

SIZE	ANKLE (CB)	CALF (CC)	FIBULAR HEAD (CD)	LENGTH (A-D)
XS	18-23 cm	32-38 cm	37-43 cm	38 cm
S	21-26 cm	37-43 cm	42-48 cm	39 cm
M	23-28 cm	42-48 cm	45-51 cm	41 cm
L	26-31 cm	47-53 cm	50-56 cm	42 cm
XL	28-33 cm	52-58 cm	55-61 cm	43 cm
XXL	31-36 cm	57-63 cm	60-66 cm	45 cm

FARROWHYBRID® SIZING CHART	
A' _____	Standard: 20-27cm
	Wide: 28-35cm

### JOBST® FARROW LINER

BNR CODE	SIZE	DESCRIPTION	QUANTITY
7666800	15-58 cm (G)	SILVER LINER AD (LEG) SIZE SMALL	
7666802	52-70 cm (G)	SILVER LINER AD (LEG) SIZE LARGE	



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

## FarrowWrap® Liner/Farrowhybrid® Liner Order Form

BNR CODE	DESCRIPTION	SIZE	LENGTH	COLOR	SIZING
7666400	FarrowHybrid® ADI Foot Compression 20-30mmHg	Small	Standard	Taupe	mens shoe size 3-6 womens shoe size 4-7
7666401	FarrowHybrid® ADI Foot Compression 20-30mmHg	Medium	Standard	Taupe	mens shoe size 7-9 womens shoe size 8-10
7666402	FarrowHybrid® ADI Foot Compression 20-30mmHg	Large	Standard	Taupe	mens shoe size 10-12 womens shoe size 11-13
7666403	FarrowHybrid® ADI Foot Compression 20-30mmHg	Extra Large	Standard	Taupe	mens shoe size 13-15 womens shoe size 14-16
7666412	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Small	Wide	Taupe	mens shoe size 3-6 womens shoe size 4-7
7666413	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Medium	Wide	Taupe	mens shoe size 7-9 wide womens shoe size 8-10 wide
7666414	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Large	Wide	Taupe	mens shoe size 10-12 wide womens shoe size 11-13 wide
7666415	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Extra Large	Wide	Taupe	mens shoe size 13-15 wide womens shoe size 14-16 wide
7666416	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Small	Wide	Black	mens shoe size 3-6 wide womens shoe size 4-7 wide
7666417	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Medium	Wide	Black	mens shoe size 7-9 wide womens shoe size 8-10 wide
7666418	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Large	Wide	Black	mens shoe size 10-12 wide womens shoe size 11-13 wide
7666419	FarrowHybrid® ADIW Foot Compression 20-30mmHg	Extra Large	Wide	Black	mens shoe size 13-15 wide womens shoe size 14-16 wide
7666404	FarrowHybridv ADII Foot Compression 30-40 mmHg	Small	Standard	Sand	mens shoe size 3-6 womens shoe size 4-7
7666405	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Medium	Standard	Sand	mens shoe size 7-9 womens shoe size 8-10
7666406	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Large	Standard	Sand	mens shoe size 10-12 womens shoe size 11-13
7666407	FarrowHybrid® ADII Foot Compression 30-40 mmHg	Extra Large	Standard	Sand	mens shoe size 13-15 womens shoe size 14-16
7666408	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Small	Wide	Sand	mens shoe size 3-6 wide womens shoe size 4-7 wide
7666409	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Medium	Wide	Sand	mens shoe size 7-9 wide womens shoe size 8-10 wide
7666410	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Large	Wide	Sand	mens shoe size 10-12 wide womens shoe size 11-13 wide
7666411	FarrowHybrid® ADIIW Foot Compression 30-40 mmHg	Extra Large	Wide	Sand	mens shoe size 13-15 wide womens shoe size 14-16 wide

**JOBST® FARROWHYBRID® FOOT COMPRESSION FARROWLINER™**

PLEASE NOTE: All liners available in Taupe unless otherwise noted in the above size description.