



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

# JOBST® FarrowWrap Upper Extremity

## Ready-To-Wear Arm FarrowWra• Custom

### Ready-To-Wear Arm

**Fabric**

LITE

**Color**

- Beige (00)
- Black (01)

**Length**

- Short (S)
- Regular (R)
- Long (L)

**Size**

- Small (2)
- Medium (3)
- Large (4)

**Side**

- Left (L)
- Right (R)

### Custom

**Style**

Custom

**Fabric**

- Classic
- LITE
- STRONG

**Color**

Tan

**Side**

- Left (L)
- Right (R)


### Trim-To-Fit

**Size**

- Medium (3)
- Extra Large (5)

**Side**

- Left (L)
- Right (R)



Circumference (c)	Left	Right
cG (axilla)*		
cF <sup>1</sup> (widest bicep)		
cF (distal bicep)		
cE (elbow crease)*		
cD <sup>1</sup> (widest forearm)		
cD (distal forearm)		
cC (wrist)*		

Lengths	Left	Right
lC-G (wrist to axilla)*		
lC-E (wrist to elbow)		

\* Circumference measurements needed for Ready to Wear garments  
 Please note that all measurements are needed for Custom garments

Comments

# JOBST® FarrowWrap Upper Extremity

## Ready-To-Wear Arm FarrowWrap • Custom

### FarrowWrap Ready-To-Wear Arm

Size	Wrist (cC)	Elbow (cE)	Axilla (cG)
S	14-18 cm	20-27 cm	22-31 cm
M	16-21 cm	25-34 cm	29-39 cm
L	19-25 cm	30-40 cm	32-45 cm
Length (C-G)	Short 40-43 cm		
	Regular 44-47 cm		
	Long 48-51 cm		

### FarrowWrap Trim-To-Fit Arm

Size	Wrist (cC)	Forearm (cD)	Elbow (cE)	Bicep (cF)	Axilla (cG)	Length (C-E)	Length (E-G)
M	15-31 cm	20-36 cm	20-38 cm	25-42 cm	25-44 cm	20-25 cm	15-23 cm
XL	15-35 cm	20-40 cm	25-42 cm	25-45 cm	30-48 cm	20-28 cm	

JOBST FARROWWRAP READY-TO-WEAR		
BSN Code	Description	Quantity
7665035	LITE Armpiece Long Small Left Tan	
7665000	LITE Armpiece Long Small Right Tan	
7665001	LITE Armpiece Long Medium Left Tan	
7665002	LITE Armpiece Long Medium Right Tan	
7665003	LITE Armpiece Long Large Left Tan	
7665004	LITE Armpiece Long Large Right Tan	
7665005	LITE Armpiece Regular Tan Small Left	
7665006	LITE Armpiece Regular Tan Small Right	
7665007	LITE Armpiece Regular Tan Medium Left	
7665008	LITE Armpiece Regular Tan Medium Right	
7665009	LITE Armpiece Regular Tan Large Left	
7665010	LITE Armpiece Regular Tan Large Right	
7665011	LITE Armpiece Short Tan Small Left	
7665012	LITE Armpiece Short Tan Small Right	
7665013	LITE Armpiece Short Tan Medium Left	
7665014	LITE Armpiece Short Tan Medium Right	
7665015	LITE Armpiece Short Tan Large Left	
7665016	LITE Armpiece Short Tan Large Right	
7665017	LITE Armpiece Long Black Small Left	
7665018	LITE Armpiece Long Black Small Right	
7665019	LITE Armpiece Long Black Medium Left	
7665020	LITE Armpiece Long Black Medium Right	
7665021	LITE Armpiece Long Black Large Left	
7665022	LITE Armpiece Long Black Large Right	

JOBST FARROWWRAP READY-TO-WEAR		
BSN Code	Description	Quantity
7665023	LITE Armpiece Regular Black Small Left	
7665024	LITE Armpiece Regular Black Small Right	
7665025	LITE Armpiece Regular Black Medium Left	
7665026	LITE Armpiece Regular Black Medium Right	
7665027	LITE Armpiece Regular Black Large Left	
7665028	LITE Armpiece Regular Black Large Right	
7665029	LITE Armpiece Short Black Small Left	
7665030	LITE Armpiece Short Black Small Right	
7665031	LITE Armpiece Short Black Medium Left	
7665032	LITE Armpiece Short Black Medium Right	
7665033	LITE Armpiece Short Black Large Left	
7665034	LITE Armpiece Short Black Large Right	

JOBST FARROWWRAP TRIM-TO-FIT		
BSN Code	Description	Quantity
7612000	LITE Armpiece (Left) Medium with 1 pair of silver arm liners	
7612002	LITE Armpiece (Right) Medium with 1 pair of silver arm liners	
7612001	LITE Armpiece (Left) Extra Large with 1 pair of silver arm liners	
7612003	LITE Armpiece (Right) Extra Large with 1 pair of silver arm liners	

JOBST FARROWWRAP CUSTOM CLASSIC		
BSN Code	Description	Quantity
7599500	Classic Arm	

JOBST FARROWWRAP CUSTOM STRONG		
BSN Code	Description	Quantity
7598700	STRONG Arm	

JOBST FARROWWRAP CUSTOM LITE		
BSN Code	Description	Quantity
7599100	LITE Custom Arm (Left or Right)	



Absolute Medical, 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: [orders@absolutemedical.com](mailto:orders@absolutemedical.com)  
 Our website: [absolutemedical.com](http://absolutemedical.com)