




Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____
 Fitter Last Name: _____ Fitter First Name: _____
 Date: _____

Classic Arm Sleeve AG1 Ready-to-Wear

ORDER SPECIFICATIONS			
<input type="checkbox"/> Quote Only	<input type="checkbox"/> Quote & Proceed	<input type="checkbox"/> Dealer Pricing	<input type="checkbox"/> MSRP
SHIPPING			
<i>Shipping rates may vary, depending on services requested and/or rates charged by carrier.</i>			
<input type="checkbox"/> FedEx® (2 day shipping) \$10.00 to business addresses; \$13.25 to residential addresses (Additional services may be available; contact JoViPak to discuss.)	<input type="checkbox"/> Check if shipping to a residence	<input type="checkbox"/> USPS Priority Mail Flat Rate® Small or Medium box For PO Boxes only; billed at current USPS rates. (If order does not fit, JoViPak will choose a different shipping method.)	



Polartec® Power Dry® Colors & Channeling	Side	Size		Length
<input type="checkbox"/> Black - Slimline Channeling	<input type="checkbox"/> Right Arm	<input type="checkbox"/> XSmall	<input type="checkbox"/> XSmall MAX	<input type="checkbox"/> XShort
<input type="checkbox"/> Royal Blue - Standard Channeling	<input type="checkbox"/> Left Arm	<input type="checkbox"/> Small	<input type="checkbox"/> Small MAX	<input type="checkbox"/> Short
<input type="checkbox"/> Royal Blue - Slimline Channeling		<input type="checkbox"/> Medium	<input type="checkbox"/> Medium MAX	<input type="checkbox"/> Regular
		<input type="checkbox"/> Large	<input type="checkbox"/> Large MAX	<input type="checkbox"/> Long
		<input type="checkbox"/> XLarge	<input type="checkbox"/> XLarge MAX	<input type="checkbox"/> XLong
<input type="checkbox"/> JoViJacket - Black Nylon & Spandex Powernet (JoViJackets are recommended as they provide the additional compression needed for maximum fit & effectiveness.)				
Comments/Quantity: _____				
Additional Items				
<input type="checkbox"/> Dycem® - donning aid				
<input type="checkbox"/> Easy Slide - donning aid				
<input type="checkbox"/> Prepaid Reduction Option				



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Classic Arm Sleeve AG1 Ready-to-Wear

All measurements are provided in centimeters.

Classic Arm Sleeve AG1 Size Chart			
	Least Wrist Circumference	Least Elbow Circumference	Axilla Circumference
	"C" Landmark	"E" Landmark	"G" Landmark
XSmall	14 - 15.5	20 - 23	23 - 27.5
XSmall MAX		22 - 25	27 - 31
Small	15.5 - 16.5	23 - 26	27.5 - 32
Small MAX		25 - 28	31 - 36
Medium	16.5 - 18	26 - 29	32 - 36
Medium MAX		28 - 31	36 - 40
Large	18 - 19	29 - 33	36 - 41
Large MAX		31 - 35	40 - 44.5
XLarge	19 - 20.5	33 - 37.5	41 - 45
XLarge MAX		35 - 39.5	44.5 - 49
<i>MAX sizes are 2cm larger at the elbow and 4cm larger at the axilla.</i>			

Classic Arm Sleeve AG1 Length Chart	
Wrist to Axilla (C to G) (measured with arm hanging down at side)	
XShort	38
Short	40.5
Regular	43
Long	45.5
XLong	48