


Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Size Charts For Basic Support Socks and Stockings

15-20, 20-30 & 30-40 mmHg



20-30 & 30-40 mmHg Basic Ribbed Compression Sock 4201 & 4202

	I	II	III	IV	V
cC	29-38 cm 11 1/2"-15"	34-43 cm 13 1/4"-17"	37-49 cm 14 1/2"-19 1/4"	41-53 cm 16 1/4"-21"	46-58 cm 18"-23"
cB	18-21 cm 7"-8 1/4"	21-24 cm 8 1/4"-9 1/2"	24-27 cm 9 1/2"-10 3/4"	27-31 cm 10 3/4"-12 1/4"	31-35 cm 12 1/4"-13 3/4"


PCSZ-01-28

length A-D

A-D Regular - 15 3/4"-18 1/4" / 40-46 cm

A-D Short - 13-15 3/4" / 34-40 cm

15-20, 20-30 & 30-40 mmHg Basic Compression Stockings & Pantyhose 4410, 4411, 4412



	I	II	III	IV
cG	44-66 cm 17 1/4"-26"	48-71 cm 18 3/4"-28"	52-76 cm 20 1/2"-29 3/4"	56-81 cm 22"-32"
cG with silicone	44-62 cm 17 1/4"-24 1/2"	48-67 cm 18 3/4"-26 1/4"	52-72 cm 20 1/2"-28 1/4"	56-77 cm 22"-30 1/4"
cC	28-38 cm 11"-15"	31-42 cm 12"-16 1/2"	34-46 cm 13 1/4"-18"	37-50 cm 14 1/2"-19 3/4"
cB	18-21.5 cm 7"-8 1/2"	21.5-25 cm 8 1/2"-9 3/4"	25-28 cm 9 3/4"-11"	28-31.5 cm 11"-12 1/2"

PCSZ-01-29

length A-d


A-D Regular - 15 3/4"-18 1/4" / 40-46 cm

A-D Short - 13-15 3/4" / 34-40 cm

length A-G

A-G / A-T Regular - 28 1/4"-32 3/4" / 72-83 cm

A-G / A-T Short - 24 3/4"-28 1/4" / 63-72 cm



15-20, 20-30 mmHg Basic Casual Compression Socks 4700, 4701

	I	II	III	IV
cC	28-38 cm 11"-15"	31-42 cm 12"-16 1/2"	34-46 cm 13 1/4"-18"	37-50 cm 14 1/2"-19 3/4"
cB	18-21.5 cm 7"-8 1/2"	21.5-25 cm 8 1/2"-9 3/4"	25-28 cm 9 3/4"-11"	28-31.5 cm 11"-12 1/2"

PCSZ-01-43

length A-d

A-D Regular - 15 3/4"-18 1/4" / 40-46 cm