



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com


Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

Size Chart For Juzo® Compression Wrap Calf Segment

Note: This size chart is not intended to replace a professional fitting by a qualified fitter.



	S	M	L	XL	2XL	3XL
Max Calf	35-41 cm 13 3/4"-16 1/2"	40-46 cm 15 3/4"-18"	44-51 cm 17 1/4"-20"	48-55 cm 19"-21 3/4"	53-60 cm 20 3/4"-23 1/2"	60-70 cm 23 1/2"-27 1/2"
Calf	29-35 cm 11 1/2"-13 3/4"	34-40 cm 13 1/4"-15 3/4"	37-44 cm 14 1/2"-17 1/4"	41-48 cm 16 1/4"-19"	46-53 cm 18"-20 3/4"	50-60 cm 19 3/4"-23 1/2"
Ankle	18-21 cm 7"-8 1/4"	21-24 cm 8 1/4"-9 1/2"	24-27 cm 9 1/2"-10 3/4"	27-31 cm 10 3/4"-12 1/4"	31-35 cm 12 1/4"-13 3/4"	36-47 cm 14"-18 1/2"

AD Regular < 40cm / 15 3/4" • AD Long > 40cm / 15 3/4"