




Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_  
 Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_  
 Date: \_\_\_\_\_

## Size Charts For Dynamic Max Compression Stockings 20-30 mmHg and 30-40 mmHg

### 20-30, 30-40 mmHg Dynamic Max Compression Stockings 3511 / 3512



**cG**  
**cC**  
**cB**

	III	IV	V
<b>cG</b>	60-82 cm 23 3/4"-32 1/4"	62-87 cm 24 1/2"-34 1/4"	77-94 cm 30 1/4"-37"
<b>cC</b>	49-55 cm 19 1/4"-21 3/4"	53-60 cm 21 1/4"-23 1/2"	55-66 cm 21 3/4"-26"
<b>cB</b>	24-27 cm 9 1/2 "-10 3/4"	27-32 cm 10 3/4"12 3/4"	32-37 cm 12 3/4"-14 3/4"

#### length A-D

**A-D Regular** - 15 3/4-18 1/4"/ 40-46 cm  
**A-D Short** 13-15 3/4"/ 34-40 cm

#### length A-G

**A-G Regular** - 28 1/4-32 3/4"/ 72-83 cm  
**A-G Short** 24 3/4-28 1/4"/ 63-72 cm