

# Juzo Silver

## Silver and Leg Health

Today's health care providers are facing escalating costs and treatment time for lower extremity patient care. In the United States alone, over 80 million people have venous disease<sup>1</sup>. These patients are at increased risk for underlying problems such as leg or foot pain, infections, swelling, skin irritations and ulcers.

The progression of venous disease states usually involves skin breakdown and slow healing times; often leading to infections and/or cellulitis. In some cases these conditions are more severe and chronic.

Many people suffering from venous disease find relief using graduated compression stockings.

Juzo Silver stockings offer the therapeutic benefits of compression, combined with the added benefits of X-Static.

## X-Static®

X-Static® fibers are made using 99% pure silver that is bonded to the textile fiber. Today X-Static® silver fibers are used in a wide range of products including wound care dressings, sports apparel, military combat uniforms, and footwear.



# Juzo Silver

## Safe and Natural

Juzo Silver compression stockings are made using X-Static. Silver is safe, natural, and for centuries has been medically recognized for its antimicrobial properties. Juzo Silver garments do not use chemicals to provide antimicrobial properties.

## Ionic Release

The antimicrobial effect of silver is generated by the release of silver ions. A release rate of at least one part per million (PPM) is required to offer effective antimicrobial efficacy. Juzo® Silver products offer a release rate well above the minimum requirement for antimicrobial activity<sup>1</sup> and solely protect the finished product from microbial growth.

## Laboratory Tested-Odor free

Juzo Silver compression stocking were tested by North American Science Associates (NAMSA) for antimicrobial activity. The tests show that Juzo Silver compression fabric protected itself from Staphylococcus Aureus in less than 4 hours<sup>4</sup>, helping to reduce odor from the garment.




Methicillin-resistant Staphylococcus Aureus (MRSA) is a type of bacteria that is resistant to certain antibiotics. These antibiotics include methicillin and other more common antibiotics such as oxacillin, penicillin and amoxicillin. Laboratory tests show that Juzo Silver products are highly effective at protecting the finished product from MRSA<sup>4</sup>.

## Wash Resistant

Juzo Silver is permanently bonded to the textile fibers. Laboratory tests show that Juzo Silver compression stockings maintain high efficacy even after 100 washes<sup>3, 4</sup>.

## Heat Dissipating

Silver is one of the most thermally conductive elements in the environment. Heat between the skin and garment is dissipated into the environment via the silver fibers. This promotes a more pleasant wearing experience.



\*Antimicrobial protects solely the finished garment from microbial growth & odor

## For Use

Juzo Silver compression stockings are recommended for the following:


- History of repeated infections or cellulitis
- Patients with a high risk of developing venous ulcers<sup>5</sup>
- Early signs of skin irritations or breakdown
- Used over top wound dressings of an active ulcer
- Patients with poor personal hygiene



## Juzo® Silver Soft

Juzo Silver Soft stockings have high elasticity and softness making them easy and comfortable to wear. Our soft compression stockings is knitted with our unique FiberSoft™ weave, and X-Static.


- Model 2060 Series
- Compression 20-30 & 30-40 mmHg
- Styles Knee-High, Thigh-High, Pantyhose
- Sizing I-V Standard, Short Length on Some Models



## Juzo® Silver Dynamic

Juzo Silver Dynamic is durable and versatile. This opaque compression stocking features a seamless, two-way stretch knit for a contoured, yet comfortable fit. Dynamic's dense knit also provides good containment. Made with X-Static®.

- Model 3560 Series
- Compression 20-30, 30-40 & 40-50 mmHg
- Styles Knee-High, Thigh-High, Pantyhose, Armsleeves & Stump Shrinkers
- Sizing Custom Made, Standard Stump Shrinkers



## Juzo® Silver Strong & Silver Expert

Juzo Silver Strong garments features X-Static's silver and a dense, micro-massaging knit that provides excellent containment. Each garment is custom made to provide an optimal fit.

- Model 3050SV (Silver Strong), 3020SV (Silver Expert )
- Compression 18-21, 23-32, 34-46, 50+ mmHg
- Styles Knee-High, Thigh-High, Pantyhose, Armsleeves
- Sizing Custom Made

# Juzo Silver

## Juzo® Silver Sole

Juzo Silver Sole Sock was designed for those with sensitive feet. A pillowed sole adds cushion to reduce blistering and callous buildup, and a channeled toe seam helps reduce toe irritation. Made with X-Static®.


- Model 5760 Series
- Compression Mild
- Styles Low Cut, Crew Length & Knee-High
- Sizing S, M, L, XL



## SilverSeal® Tubular Component

SilverSeal X-Static® Tubular Component is used to provide a non-compressive protective layer between the skin and short stretch bandages. Made with X-Static®, this tubular component offers a patient the added benefits of silver early during therapy.

- Model 1600 Series
- Styles Cut to Fit Upper or Lower Extremities
- Sizing 2" x 25 yards, 3" x 25 yards, 4" x 25 yards and 6" x 25 yards



## Silver Ulcer Liner

The Juzo Silver Ulcer Liner, made with X-Static silver, offers mild compression and can be worn by itself or used as an under layer for a Juzo compression stocking

- Model 9641 Series
- Compression Mild
- Styles Knee High
- Sizing S, L



### Foot Notes

- 1 American College of Phlebology
  - 2 American Diabetes Association
  - 3 Testing performed by Noble Fiber Technologies Advanced Material Development Laboratory. Data on file.
  - 4 Antimicrobial testing performed by North American Science Associates. Data on file.
  - 5 Samson RH, Showalter DP, Stockings and the prevention of recurrent venous ulcers, Dermatologic Surgery, 1996 Apr; 22(4):373-6
  - 6 Department of Health & Human Services, CMS, Program Memorandum AB-03-090, June 20,2003
- © X-Static and X-Static - The Silver Fiber are registered trademarks of Noble Fiber Technologies.  
© Juzo is a registered trademark of Julius Zorn, Inc.  
© 2016 - Julius Zorn, Inc

Juzo U.S.A.  
www.juzousa.com  
1 800 222 4999

Juzo Canada Ltd.  
www.juzo.ca  
1 844 229 4999

