

Patient Name: _____

Fitter Name: _____

Date: _____




Absolute Medical. 1843 W Hubbard St #2A
Chicago, IL 60622
Tel: 312 233 2207 Fax: 866 860 9358
To Order Online: orders@absolutemedical.com
Our website: absolutemedical.com

Size Charts For Compression Arm Sleeves & Gauntlets

15-20mmHg, 20-30 mmHg & 30-40 mmHg

2000, 2001, 2002, 3511, 3512
20-30 and 30-40 mmHg




C-G/C-H Regular - C-G length <17" / 43 cm
C-G/C-H Long - C-G length >17" / 43 cm

	I	II	III	IV	V	VI
g max	26.5-33 cm 10 1/2"-13"	28.5-35.0 cm 11 1/4"-13 3/4"	30.5-37.5 cm 12"-14 3/4"	32.5-39.5 cm 12 3/4"-15 1/2"	34.0-42.0 cm 13 1/2"-16 1/2"	36.0-44.0 cm 14 1/2"-17 1/2"
g	21.5-26.5 cm 8 1/2"-10 1/2"	23.5-28.5 cm 9 1/4"-11 1/4"	25.0-31.0 cm 9 3/4"-12 1/4"	27.0-33.0 cm 10 1/2"-13"	29.0-35.0 cm 11 1/2"-13 3/4"	30.5-37.5 cm 12"-14 3/4"
e max	22.5-27.5 cm 8 3/4"-10 3/4"	24.5-29.5 cm 9 3/4"-11 1/2"	26-32 cm 10 1/4"-12 1/2"	27.5-33.5 cm 10 3/4"-13 1/4"	29.0-35.0 cm 11 1/2"-13 3/4"	30.0-37.0 cm 11 3/4"-14 1/2"
e	20.5-25.5 cm 8"-10"	22-27 cm 8 3/4"-10 3/4"	23.5-28.5 cm 9 1/4"-11 1/4"	25.0-30.0 cm 9 3/4"-11 3/4"	26.0-32.0 cm 10 1/4"-12 1/2"	27.5-33.5 cm 10 3/4"-13 1/4"
c	14.0-15.5 cm 5 1/2"-6"	15.5-17.0 cm 6"-6 3/4"	17.0-18.5 cm 6 3/4"-7 1/4"	18.5-20 cm 7 1/4"-7 3/4"	20.0-21.5 cm 7 3/4"-8 1/2"	21.5-23 cm 8 1/2"-9"

PCSZ-01-04b

Size Chart for Juzo® 2001 / 2002, 2301 / 2302 Gloves & Gauntlets



	1-S	2-M	3-L	4-XL
a	17.5-19.5 cm 6 7/8"-7 5/8"	19.5-21.0 cm 7 5/8"-8 1/4"	21.0-23.0 cm 8 1/4"-9 1/8"	23.0-25.5 cm 9 1/8"-10 1/8"
c	14.5-16.5 cm 5 5/8"-6 1/2"	16.5-18.0 cm 6 1/2"-7"	18.0-20.0 cm 7"-7 7/8"	20.0-22 cm 7 7/8"-8 5/8"

Size Chart for Juzo® 3021 / 3022 Gauntlet

	1		2		3		4		5		6	
	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm	X cm	Z cm
①	6.5	5	6.5	5	6.5	5	7.5	6	7.5	6	7.5	6
②	6.2	5	6.7	5.5	7.1	5.9	7.6	6.4	8	6.8	8.5	7.3
③	6	5	6.3	5.3	6.7	5.7	7.1	6	7.5	6.4	7.6	6.7
④	5.7	4	6	4.2	6.3	4.5	6.7	4.8	7.1	5.1	7.4	5.4
⑤	5.5	4	5.9	4.2	6.3	4.5	6.6	4.8	7	5.1	7.3	5.4
a	17.5 cm 6 7/8"		19 cm 7 1/2"		20.5 cm 8"		22 cm 8 3/4"		23.5 cm 9 1/4"		25 cm 9 3/4"	
c	14.5 cm 5 3/4"		16 cm 6 1/4"		17.5 cm 7"		19 cm 7 1/2"		20.5 cm 8"		22 cm 8 3/4"	

PCSZ-01-06b