



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

### circaid capri

<i>size</i>	<i>circumference in inches</i>		<i>circumference in centimeters</i>	
	<i>waist (T)</i>	<i>thigh (G)</i>	<i>waist (T)</i>	<i>thigh (G)</i>
<b>A</b>	25 - 31	18 - 24	63 - 79	46 - 62
<b>B</b>	26 - 33	19 - 26	67 - 85	49 - 65
<b>C</b>	28 - 36	20 - 27	72 - 91	52 - 70
<b>D</b>	30 - 38	22 - 29	76 - 96	55 - 73
<b>E</b>	32 - 40	24 - 32	81 - 102	60 - 81
<b>F</b>	33 - 43	25 - 33	85 - 108	63 - 84
<b>G</b>	35 - 44	26 - 35	90 - 113	67 - 90
<b>H</b>	37 - 47	27 - 37	94 - 119	70 - 93
<b>I</b>	39 - 49	28 - 38	99 - 125	72 - 96
<b>J</b>	41 - 51	30 - 40	103 - 130	77 - 102

When taking measurements, measure at the largest circumference point. When selecting a size, choose the smallest size match for the most snug fit. Choose the largest size match for a more relaxed, easy-to-wear fit.