




Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel: 312 233 2207 Fax: 866 860 9358  
 To Order Online: orders@absolutemedical.com  
 Our website: absolutemedical.com

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_



**upper leg kit:**

\*cG cannot be more than 30cm larger than cE2.

length			Width	
short	Standard	Long	Regular	Wide
D1-G	D1-G	D1-G	cG < 80cm	cG < 125cm
< 37cm	37 - 41cm	>41cm	cE2 < 65cm	cE2 < 110cm

**knee kit: (the reduction knee kit is only available in one length and width)**

\*cE2 cannot be more than 30cm larger than cC.

length	Width
30cm	cE 2 < 110cm cE < 102cm cC < 80cm

**lower leg kit:**

\*cC cannot be more than 30cm larger than cB.

length			Width	
short	Standard	Long	Regular	Wide
l mall - D1	l mall - D1	l mall - D1	cC < 65cm	cC < 80cm
< 32cm	32 - 46cm	>36cm	cB < 55cm	cB < 65cm

**arm kit:**

\*If cB is >25cm hand wrap will not fit through reduction arm component may still be used.

length		Width	
standard	Long	Regular	Wide
l C-G < 45cm	l C-G > 45cm	cF < 50cm cB < 25cm*	cF 50-80xm cB < 25cm*

