



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

mediven comfort 15-20, 20-30, 30-40 mmHg size chart

<i>circumference</i>	I	II	III	IV	V	VI	VII
panty & thigh	16½" – 22½" 42–57 cm	17¾" – 24¼" 45–61.5 cm	19¼" – 26¼" 49–66.5 cm	20½" – 28¾" 52–73 cm	22" – 31" 56–79 cm	23¾" – 33¼" 60–84.5 cm	25½" – 34½" 64–87.5 cm
extra-wide calf	N/A	N/A	16½" – 19½" 42–49.5 cm	17¾" – 20¾" 45–53 cm	19" – 22" 48 – 56 cm	20¼" – 23¼" 51.5–59 cm	21½" – 24½" 54.5– 62 cm
calf	11" – 14" 28–35.5 cm	11¾" – 15¼" 30–39 cm	12¾" – 16½" 32.5–42 cm	13¾" – 17¾" 35–45 cm	14½" – 19" 37–48.5 cm	15½" – 20¼" 39.5–51.5 cm	16¼" – 21½" 41.5–54.5 cm
ankle	7" – 7¾" 18–20 cm	8" – 8¾" 20.5–22.5 cm	9" – 9¾" 23–25 cm	10" – 10¾" 25.5–27.5 cm	11" – 11¾" 28–30 cm	12" – 13" 30.5–33 cm	13¼" – 14¼" 33.5–36 cm
leg length	calf ▶ standard: 15"–19" (38–47.5 cm) petite: 12–15" (31–38 cm) panty & thigh ▶ standard: 28"–34" (71–86.5 cm) petite: 23"–28" (58.5–71 cm)						