



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

mediven harmony sleeves, sleeve with gauntlet, gauntlet and gloves*

| <i>circumference (cm)</i> | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
|---------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| "G" extra-wide | 28–31 | 30–34 | 32.5–36.5 | 34.5–39.5 | 36.5–41.5 | 39–44 | 41–46.5 | 44.5–50.5 |
| "G" standard | 25–28 | 27.5–31 | 29.5–33.5 | 31.5–36 | 33–37.5 | 35.5–40.5 | 37.5–42.5 | 40–45.5 |
| "E" extra-wide | 24–26.5 | 26–29.5 | 28–31.5 | 29.5–34 | 31–35.5 | 33–37.5 | 35–40 | 38–43 |
| "E" standard | 23–25 | 24–27 | 25.5–29 | 27–31 | 28.5–32.5 | 30–34.5 | 32–36.5 | 34.5–39.5 |
| "C1" | 15–16.5 | 16–18.5 | 17.5–20 | 19–22 | 20.5–24 | 22.5–26 | 24.5–28.5 | 27–31 |
| "C" | 13.5–14.5 | 14–16 | 15.5–18 | 17–20 | 18.5–21.5 | 20–23.5 | 22–25.5 | 24.5–28.5 |
| "A" | 15–18 | 15.5–19.5 | 16.5–20.5 | 18.5–23 | 19.5–24 | 21.5–26 | 23.5–28 | 25.5–30 |

gloves & gauntlets available in sizes 2-7

armsleeves available in sizes 1-8

*All measurements in centimeters.