



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

mediven sheer & soft 15–20, 20–30, 30–40 mmHg size chart

<i>circumference</i>	I	II	III	IV	V	VI	VII
panty & thigh	16½"– 22½" 42–57 cm	17¾"– 24¼" 45–61.5 cm	19¼"– 26¼" 49–66.5 cm	20½"– 28¾" 52–73 cm	22"– 31" 56–79 cm	23¾"– 33¼" 60–84.5 cm	25½"– 34½" 64–87.5 cm
calf	11"– 14" 28–35.5 cm	11¾"– 15¼" 30–39 cm	12¾"– 16½" 32.5–42 cm	13¾"– 17¾" 35–45 cm	14½"– 19" 37–48.5 cm	15½"– 20¼" 39.5–51.5 cm	16¼"– 21½" 41.5–54.5 cm
ankle	7"– 7¾" 18–20 cm	8"– 8¾" 20.5–22.5 cm	9"– 9¾" 23–25 cm	10"– 10¾" 25.5–27.5 cm	11"– 11¾" 28–30 cm	12"– 13" 30.5–33 cm	13¼"– 14¼" 33.5–36 cm
leg length	calf ▶ standard: 15"–19" (38–47.5 cm) petite: 12–15" (31–38 cm) panty & thigh ▶ standard: 28"–34" (71–86.5 cm) petite: 23"–28" (58.5–71 cm)						