



Absolute Medical. 1843 W Hubbard St #2A  
Chicago, IL 60622  
Tel: 312 233 2207 Fax: 866 860 9358  
To Order Online: [orders@absolutemedical.com](mailto:orders@absolutemedical.com)  
Our website: [absolutemedical.com](http://absolutemedical.com)

Patient Last Name: \_\_\_\_\_ Patient First Name: \_\_\_\_\_

Fitter Last Name: \_\_\_\_\_ Fitter First Name: \_\_\_\_\_

Date: \_\_\_\_\_

## CHIPSLEEVE ARM

### PRODUCT INFORMATION


Arm \_\_\_\_\_ Size \_\_\_\_\_ Length \_\_\_\_\_ Item # \_\_\_\_\_

☐ LEFT ☐ RIGHT

OVERSLEEVE Color

☐ BLACK ☐ NAVY ☐ PINK ☐ TIE DYE ☐ BLUE HORIZON ☐ LEOPARD

**NOTE:** If no color is selected, BLACK will be included.



### CHIPSLEEVE ARM SIZING CHART (A-C)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C	23-32cm	28-37cm	33-43cm	39-49cm	44-54cm
B1	22-30cm	26-34cm	30-39cm	35-44cm	39-49cm
B	20-27cm	24-31cm	28-35cm	32-39cm	36-43cm
A	14-17cm	15-18cm	16-19cm	18-20cm	19-22cm

#### LEFT

SHORT	2631-ARS-L	2632-ARS-L	2633-ARS-L	2634-ARS-L	2635-ARS-L
REGULAR	2631-ARR-L	2632-ARR-L	2633-ARR-L	2634-ARR-L	2635-ARR-L
LONG	2631-ARL-L	2632-ARL-L	2633-ARL-L	2634-ARL-L	2635-ARL-L

#### RIGHT

SHORT	2631-ARS-R	2632-ARS-R	2633-ARS-R	2634-ARS-R	2635-ARS-R
REGULAR	2631-ARR-R	2632-ARR-R	2633-ARR-R	2634-ARR-R	2635-ARR-R
LONG	2631-ARL-R	2632-ARL-R	2633-ARL-R	2634-ARL-R	2635-ARL-R

Custom sizes and configurations are available for CHIPSLEEVE.  
For more information, please speak with your Absolute  
Sales Representative or Customer Care: 312-233-2207