



Absolute Medical. 1843 W Hubbard St #2A  
 Chicago, IL 60622  
 Tel. 312 233 2207 Fax 866 860 9358  
 To Order Online: orders@absolutemedical.net  
 Our website: absolutemedical.com

1BØBT\BNF \_\_\_\_\_  
 ØSBT\BNF \_\_\_\_\_  
 ØEJ \_\_\_\_\_

ØS First Name: \_\_\_\_\_

# COMPREFIT PLUS

## PRODUCT INFORMATION

Leg	Size	Length	Item #	CompreBoot Size
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				
Leg	Size	Length	Item #	CompreBoot Size
<input type="checkbox"/> LEFT <input type="checkbox"/> RIGHT				

### COMPREFIT PLUS REGULAR SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
C	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
B	19-29cm	24-34cm	28-38cm	36-46cm	39-50cm
A	15-24cm	16-25cm	20-31cm	27-37cm	28-38cm
BLACK					
REGULAR	1101-BPR	1102-BPR	1103-BPR	1104-BPR	1105-BPR

### COMPREFIT PLUS TALL SIZING CHART (A-C1)

	SMALL	MEDIUM	LARGE	XLARGE	XXLARGE
C1	24-34cm	29-39cm	34-44cm	43-53cm	50-60cm
C	21-31cm	26-36cm	30-40cm	39-49cm	45-55cm
B	16-26cm	20-30cm	25-35cm	31-41cm	35-45cm
A	15-24cm	16-25cm	20-31cm	27-37cm	28-38cm
BLACK					
TALL	1101-BPT	1102-BPT	1103-BPT	1104-BPT	1105-BPT

### COMPREBOOT SIZING CHART (I-K)

	SMALL		MED / LARGE		XL / XXL	
	REGULAR	LONG	REGULAR	LONG	REGULAR	LONG
I	28-36cm	28-36cm	≤39cm	≤39cm	≤44cm	≤44cm
J	22-26cm	22-26cm	≤30cm	≤30cm	≤33cm	≤33cm
K	14-18cm	19-23cm	18-20cm	23-25cm	20-22cm	25-27cm

### CIRCUMFERENCE

LEFT RIGHT

\_\_\_\_\_ C1

\_\_\_\_\_ C

\_\_\_\_\_ B

\_\_\_\_\_ A

\_\_\_\_\_ Ø POINT

\_\_\_\_\_ I

\_\_\_\_\_ J

\_\_\_\_\_ K: Measure **length** on medial side of foot from heel to 1st metatarsal head

