



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel: 312 233 2207 Fax: 866 860 9358
 To Order Online: orders@absolutemedical.com
 Our website: absolutemedical.com

Patient Last Name: _____ Patient First Name: _____

Fitter Last Name: _____ Fitter First Name: _____

Date: _____

HAND COMPRESSION WEAR

PRODUCT INFORMATION


Hand Style

<input type="checkbox"/> MEDAHAND	Size	Item #
<input type="checkbox"/> MEDAGLOVE	Size	Item #
<input type="checkbox"/> DORSAL POCKET GLOVE	Size	Item #

CIRCUMFERENCES

Measure circumference at Palm and record on line E.

Measure at Wrist, and record on line A.



HAND SIZING CHART (A-E)

	SMALL	MEDIUM	LARGE	X LARGE
E	15-19cm	19-22cm	22-26cm	26-29cm
A	14-20cm	16-22cm	18-24cm	20-26cm

MEDAHAND

RIGHT	1201-HDR-R	1202-HDR-R	1203-HDR-R	1204-HDR-R
LEFT	1201-HDR-L	1202-HDR-L	1203-HDR-L	1204-HDR-L

MEDAGLOVE

RIGHT	1201-GLR-R	1202-GLR-R	1203-GLR-R	1204-GLR-R
LEFT	1201-GLR-L	1202-GLR-L	1203-GLR-L	1204-GLR-L

DORSAL POCKET GLOVE

RIGHT	1201-DGR-R	1202-DGR-R	1203-DGR-R	1204-DGR-R
LEFT	1201-DGR-L	1202-DGR-L	1203-DGR-L	1204-DGR-L