



Absolute Medical. 1843 W Hubbard St #2A
 Chicago, IL 60622
 Tel. 312 233 2207 Fax 866 860 9358
 To Order Online: orders@absolutemedical.net
 Our website: absolutemedical.com

1 ~~BO~~ ~~NT~~ ~~BNF~~ _____
~~ST~~ ~~BNF~~ _____ Fitter First Name: _____
~~ST~~ _____


NECK/MANDIBLE PAD

PRODUCT INFORMATION


Size

Item #

SMALL MEDIUM LARGE X-LARGE



HEAD - measure circumference starting at the top of the head down to under the chin, back to top of head.



NECK - measure circumference at largest part of neck.

NECK/MANDIBLE PAD SIZING CHART

	SMALL	MEDIUM	LARGE	XLARGE
HEAD	46-56cm	52-62cm	56-66cm	58-68cm
NECK	28-33cm	33-38cm	38-43cm	43-58cm
	2231-NM	2232-NM	2233-NM	2234-NM