




WASH AND CARE GUIDE FOR COMPRESSION GARMENTS

In order to prevent damage to the material of the compression garment, avoid sticking fingernails into the knitted material and do not wear jewelry while putting on or taking off the garment.

FLAT KNIT



-  Machine Washable
-  Mild Detergent
-  Use Lukewarm Water
-  Tumble Dry Low
-  **NO Fabric Softener or Woolite**
-  **DO NOT Bleach**
-  **DO NOT Iron**

READY-TO-WEAR



-  Machine Washable
-  Mild Detergent
-  Use Lukewarm Water
-  Tumble Dry Low
-  **NO Fabric Softener or Woolite**
-  **DO NOT Bleach**
-  **DO NOT Iron**

VELCRO



-  Hand Wash
-  Mild Detergent
-  Use Lukewarm Water
-  Hang or Lay Flat
-  **DO NOT Submerge in Water**
-  **DO NOT Bleach**
-  **DO NOT Iron**

NIGHTTIME



-  Machine Washable
-  Mild Detergent
-  Use Lukewarm Water
-  Tumble Dry Low
-  **NO Fabric Softener**
-  **DO NOT Bleach**
-  **DO NOT Iron**

BANDAGES



-  Hand Wash
-  Mild Detergent
-  Use Lukewarm Water
-  Dry Flat
-  **DO NOT Machine Dry**
-  **NO Detergent**
-  **DO NOT Iron**

PUMP SLEEVES



-  Hand Wash
-  Mild Detergent
*Applied with soft cloth
-  Air Dry or Lay Flat
-  **DO NOT Use Scrubbing Pads**
-  **DO NOT Bleach**
-  **DO NOT Submerge in Water**

Follow Us! @absolutemedical

