


MOBIDERM Technology

MOBIDERM is a medical device comprised of foam cubes encased between two non-woven bandages.

MOBIDERM Technology can be used under a reducing bandage (in its bandage or pad forms), or incorporated into mobilizing garments.

The pressure difference created between the support area of the blocks and the surrounding area stimulates the flow of lymphatic fluid, optimizing drainage efficacy.⁽³⁾

PATENTED⁽⁴⁾



Two sizes of foam blocks adapt to the various needs of the affected area



**SMALL
BLOCKS**

- More effective in highly induced areas
- Better for small morphologies or body areas: pediatric patients, hands/foot, breast, and more.



**BIG
BLOCKS**

- Better adapted along the limb

MOBIDERM autofit garments are recommend for patients with all stages of lymphedema, dysmorphic or regular morphology, progressive edema.

MOBIDERM Standard

MOBIDERM garments are also available in non adjustable version.



Mitten



Glove




Armsleeve



Phone: (312) 233-2207
Fax: (866) 860-9358
www.absolutemedical.com


ORDERING INFORMATION

MOBIDERM autofit
Below-Knee




Size	Ankle Circumference (at the narrowest point) (cB)	Calf Circumference (at the widest point) (cC)	Below-Knee Circumference (on the head of the fibula) (cD)	Short ℓD = 38-40 cm	Normal ℓD = 40-42 cm	Long ℓD = 42-44 cm	Xtra Long ℓD = 44-46 cm
1	16-24 cm	24-35 cm	25-36 cm	373402 20100231	373402 20100232	373402 20100233	373402 20100234
2	20-28 cm	32-43 cm	32-43 cm	373402 20200231	373402 20200232	373402 20200233	373402 20200234
3	24-32 cm	40-51 cm	39-50 cm	373402 20300231	373402 20300232	373402 20300233	373402 20300234
4	28-36 cm	48-59 cm	46-57 cm	373402 20400231	373402 20400232	373402 20400233	373402 20400234

MOBIDERM autofit
Thigh-High




Size	Ankle Circumference (at the narrowest point) (cB)	Calf Circumference (at the widest point) (cC)	Thigh Circumference (at the widest point) (cG)	Short ℓG = 67-72 cm	Normal ℓG = 72-77 cm	Long ℓG = 77-82 cm	Xtra Long ℓG = 82-87 cm
1	16-24 cm	24-35 cm	46-57 cm	373502 20100231	373502 20100232	373502 20100233	373502 20100234
2	20-28 cm	32-43 cm	53-65 cm	373502 20200231	373502 20200232	373502 20200233	373502 20200234
3	24-32 cm	40-51 cm	61-72 cm	373502 20300231	373502 20300232	373502 20300233	373502 20300234
4	28-36 cm	48-59 cm	68-80 cm	373502 20400231	373502 20400232	373502 20400233	373502 20400234

MOBIDERM autofit
Armsleeve




Size	Hand Circumference (cA)	Wrist Circumference (cC)	Forearm Circumference (at the widest point) (cD)	Axillary Arm Circumference (cG)	Normal ℓ = 40-45 cm	Long ℓ = 45-50 cm
1	17-21 cm	15-19 cm	22-28 cm	25-32 cm	373303 20100261	373303 20100271
2	17-21 cm	16-20 cm	24-30 cm	29-36 cm	373303 20200261	373303 20200271
3	18-22 cm	17-21 cm	26-32 cm	33-40 cm	373303 20300261	373303 20300271
4	18-22 cm	17-21 cm	28-34 cm	37-44 cm	373303 20400261	373303 20400271
5	19-23 cm	18-22 cm	30-36 cm	41-48 cm	373303 20500261	373303 20500271
6	20-24 cm	19-23 cm	31-37 cm	45-52 cm	373303 20600261	373303 20600271

MOBIDERM Standard
Armsleeve



Size	Wrist Circumference (cC)	Elbow Circumference (cE)	Axillary Arm Circumference (cG)	Normal ℓ = 40-45 cm	Long ℓ = 45-50 cm
1	14-16 cm	22-25 cm	25-29 cm	373003 20100232	373003 20100233
1+	14-16 cm	24-27 cm	28-32 cm	373003 40100232	373003 40100233
2	16-18 cm	25-28 cm	29-33 cm	373003 20200232	373003 20200233
2+	16-18 cm	27-30 cm	31-35 cm	373003 40200232	373003 40200233
3	18-20 cm	28-31 cm	32-36 cm	373003 20300232	373003 20300233
3+	18-20 cm	30-33 cm	35-39 cm	373003 40300232	373003 40300233
4	20-22 cm	31-34 cm	35-39 cm	373003 20400232	373003 20400233
4+	20-22 cm	33-36 cm	38-42 cm	373003 40400232	373003 40400233
5	22-24 cm	34-37 cm	38-42 cm	373003 20500232	373003 20500233
5+	22-24 cm	36-39 cm	41-45 cm	373003 40500232	373003 40500233

MOBIDERM Standard
Glove & Mitten



Size	Wrist Circumference (cC)	Hand Circumference (cA)	Glove		Mitten	
			Right	Left	Right/Left	Right/Left
1	14-16 cm	15-17 cm	373201 20100201	373201 20100202	373101 20100299	
2	16-18 cm	17-19 cm	373201 20200201	373201 20200202	373101 20200299	
3	17-19 cm	19-21 cm	373201 20300201	373201 20300202	373101 20300299	
4	19-21 cm	21-23 cm	373201 20400201	373201 20400202	373101 20400299	
5	20-22 cm	23-25 cm	373201 20500201	373201 20500202	373101 20500299	
6	21-23 cm	25-27 cm	373201 20600201	373201 20600202	373101 20600299	

(1) Mestre - Interest of an auto-adjustable night-time compression sleeve (MOBIDERM Autofit) in maintenance phase of upper limb lymphoedema: the MARILYN pilot RCT - Support Care Cancer (25:2455-2462) - 2017 - Page 7

(2) Quéré - Prospective multicenter observational study of lymphoedema therapy: POLIT study - Vascular Disease Journal - 2014 - Pages 4 and 5

(3) Todd - Mobiderm Autofit: an adjustable sleeve that enables patients to self-manage lymphoedema - Chronic oedema - April 2018 - Page 6

(4) Only in following countries: Algeria, Austria, Belgium, Canada, China, Czech Republic, Denmark, Finland, France, Germany, Hungary, India, Morocco, Italy, Netherlands, Poland, Portugal, Romania, Russian Federation, Slovakia, Slovenia, Spain, Sweden, Switzerland, Tunisia, United-Kingdom

Biflex - Indications: Venous insufficiency, simple and complicated varicose veins, reduction of venous and posttraumatic edema, dermatitis/ lipodermatosclerosis, varicose ulcers, deep vein thrombosis (prevention and confirmed accident), post-thrombotic syndrome, following sclerotherapy, treatment of lymphedema. Contraindications: PAD (peripheral arterial disease), skin wounds, abscess, infected skin conditions, extra-anatomic bypass, uncontrolled heart failure.

Biflexideal® - Indications: Treatment of lymphedema, simple and complicated varicose veins, following sclerotherapy, following surgery, reduction of post-traumatic edema, varicose ulcers. Contraindications: Do not apply directly to an open wound.

Mobiderm Autofit and Standard - Indications: Maintenance of volume reduction in lymphedema. Contraindications: Skin infection of the limb or acute inflammation, known allergy to components used, septic thrombosis, severe peripheral neuropathy of the limb. Lower limb - PAD (peripheral arterial disease) of lower limbs with ABI < 0.6, congestive heart failure, advanced diabetic microangiopathy, phlegmasia cerulea dolens, extra-anatomic bypass. Upper limb - pathology of brachial plexus, vasculitis of the extremities.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. All the medical devices mentioned on this document are CE marked according to the European council directive 93/42/EEC and its relatives or the Regulation 2017/745 on medical devices, unless specifically identified as "NOT CE MARKED". Products mentioned in this document are CE class I devices. Please contact Thuasne should you need any additional information on devices classification.

Please read carefully the instructions for use, indications and contraindications of the product. Last revision date: 2020/07. Ref. 002



- Order
- Quotation
- Renewal

Patient's First Name: _____

Patient's Last Name: _____

Gender: M F Pediatric

Patient Height: _____

Case No. for renewal: _____

1st Treatment

Left Leg Right Leg

Fill out one for each side

MOBIDERM Lower Extremity

Auto Fit

Comments: _____

Date: _____ Quantity: _____

