

ExoStrong Thigh High 54-107

All sizes in centimeters

Short	Foot Length: up to 17		Thigh Length: up to 71	
Description	Ankle Girth	Calf Girth	Thigh Girth	Item Code
Small	18 - 20	27 - 36	41 - 55	54-107-44-XX-3220
Medium	20 - 24	29 - 39	45 - 61	54-107-45-XX-3220
Large	24 - 28	32 - 43	50 - 66	54-107-46-XX-3220
X-Large	28 - 31	33 - 45	56 - 71	54-107-47-XX-3220
XX-Large	32 - 34	37 - 50	66 - 81	54-107-48-XX-3220
Average	Foot Length: up to 19		Thigh Length: up to 74	
Description	Ankle Girth	Calf Girth	Thigh Girth	Item Code
Small	18 - 20	27 - 36	41 - 55	54-107-54-XX-3220
Medium	20 - 24	29 - 39	45 - 61	54-107-55-XX-3220
Large	24 - 28	32 - 43	50 - 66	54-107-56-XX-3220
X-Large	28 - 31	33 - 45	56 - 71	54-107-57-XX-3220
XX-Large	32 - 34	37 - 50	66 - 81	54-107-58-XX-3220
Tall	Foot Length: up to 22		Thigh Length: up to 77	
Description	Ankle Girth	Calf Girth	Thigh Girth	Item Code
Small	18 - 20	27 - 36	41 - 55	54-107-64-XX-3220
Medium	20 - 24	29 - 39	45 - 61	54-107-65-XX-3220
Large	24 - 28	32 - 43	50 - 66	54-107-66-XX-3220
X-Large	28 - 31	33 - 45	56 - 71	54-107-67-XX-3220
XX-Large	32 - 34	37 - 50	66 - 81	54-107-68-XX-3220

ExoStrong Below Knee 54-104

All sizes in centimeters

Short	Foot Length: up to 17		Calf Length: up to 38	
Description	Ankle Girth	Calf Girth	Calf Girth	Item Code
Small	18 - 20		27 - 36	54-104-44-XX-32X0
Medium	20 - 24		29 - 39	54-104-45-XX-32X0
Large	24 - 28		32 - 43	54-104-46-XX-32X0
X-Large	28 - 31		33 - 45	54-104-47-XX-32X0
XX-Large	32 - 34		37 - 50	54-104-48-XX-32X0
Average	Foot Length: up to 19		Calf Length: up to 41	
Description	Ankle Girth	Calf Girth	Calf Girth	Item Code
Small	18 - 20		27 - 36	54-104-54-XX-32X0
Medium	20 - 24		29 - 39	54-104-55-XX-32X0
Large	24 - 28		32 - 43	54-104-56-XX-32X0
X-Large	28 - 31		33 - 45	54-104-57-XX-32X0
XX-Large	32 - 34		37 - 50	54-104-58-XX-32X0
Tall	Foot Length: up to 22		Calf Length: up to 44	
Description	Ankle Girth	Calf Girth	Calf Girth	Item Code
Small	18 - 20		27 - 36	54-104-64-XX-32X0
Medium	20 - 24		29 - 39	54-104-65-XX-32X0
Large	24 - 28		32 - 43	54-104-66-XX-32X0
X-Large	28 - 31		33 - 45	54-104-67-XX-32X0
XX-Large	32 - 34		37 - 50	54-104-68-XX-32X0