



# circaid<sup>®</sup> profile

## ready-to-wear options

- Energy oversleeve is available in four color choices.
- Offers the widest ready-to-wear size ranges in the market with 16 sizes and 2 lengths.
- For easy ordering, uses same size chart as mediven<sup>®</sup> harmony and comfort arm sleeves

### Available Colors

| Item   | Description                                       | MIDNIGHT | GREY | MAGENTA | BLUE GIRAFFE | Size   |
|--------|---|----------|------|---------|--------------|--------|
|        |   | 5        | 2    | M       | 3            |        |
| CFAK11 | Arm sleeve with Energy Oversleeve                 | •        | •    | •       | •            | II-IX  |
| CFAK21 | Arm sleeve with Energy Oversleeve Extra Wide      | •        | •    | •       | •            | III-IX |
| CFAK12 | Arm sleeve with Energy Oversleeve Long            | •        | •    | •       | •            | II-IX  |
| CFAK22 | Arm sleeve with Energy Oversleeve Long Extra Wide | •        | •    | •       | •            | III-IX |

### Available Colors

| Item   | Description                       | MIDNIGHT | GREY | MAGENTA | BLUE GIRAFFE | Size   |
|--------|-----------------------------------|----------|------|---------|--------------|--------|
|        |                                   | 5        | 2    | M       | 3            |        |
| CFAS11 | Energy Oversleeve                 | •        | •    | •       | •            | II-IX  |
| CFAS21 | Energy Oversleeve Extra Wide      | •        | •    | •       | •            | III-IX |
| CFAS12 | Energy Oversleeve Long            | •        | •    | •       | •            | II-IX  |
| CFAS22 | Energy Oversleeve Extra Wide Long | •        | •    | •       | •            | III-IX |

### Standard circaid profile measuring chart

| circum. (cm)   | 2       | 3         | 4       | 5         | 6         | 7         | 8         | 9     |
|--|---------|-----------|---------|-----------|-----------|-----------|-----------|-------|
| "G" standard   | 27.5–31 | 29.5–33.5 | 31.5–36 | 33–37.5   | 35.5–40.5 | 37.5–42.5 | 40–45.5   | 43–49 |
| "E" standard   | 24–27   | 25.5–29   | 27–31   | 28.5–32.5 | 30–34.5   | 32–36.5   | 34.5–39.5 | 38–43 |
| "C1"   | 16–18.5 | 17.5–20   | 19–22   | 20.5–24   | 22.5–26   | 24.5–28.5 | 27–31     | 28–32 |
| "C"  | 14–16   | 15.5–18   | 17–20   | 18.5–21.5 | 20–23.5   | 22–25.5   | 24.5–28.5 | 25–29 |
| Arm length/Longitud del brazo/Longueur de bras: (ℓC-G:) Regular: <48 cm Long: >48 cm |         |           |         |           |           |           |           |       |

### Extra Wide circaid profile measuring chart

| circum. (cm)   | 2   | 3         | 4         | 5         | 6       | 7         | 8         | 9     |
|--|-----|-----------|-----------|-----------|---------|-----------|-----------|-------|
| "G" extra-wide   | N/A | 32.5–36.5 | 34.5–39.5 | 36.5–41.5 | 39–44   | 41–46.5   | 44.5–50.5 | 48–54 |
| "E" extra-wide   | N/A | 28–31.5   | 29.5–34   | 31–35.5   | 33–37.5 | 35–40     | 38–43     | 41–46 |
| "C1"   | N/A | 17.5–20   | 19–22     | 20.5–24   | 22.5–26 | 24.5–28.5 | 27–31     | 28–32 |
| "C"  | N/A | 15.5–18   | 17–20     | 18.5–21.5 | 20–23.5 | 22–25.5   | 24.5–28.5 | 25–29 |
| Arm length/Longitud del brazo/Longueur de bras: (ℓC-G:) Regular: <48 cm Long: >48 cm |     |           |           |           |         |           |           |       |

