

Standard circaid profile measuring chart

circum. (cm)	2	3	4	5	6	7	8	9
"G" standard	27.5–31	29.5–33.5	31.5–36	33–37.5	35.5–40.5	37.5–42.5	40–45.5	43–49
"E" standard	24–27	25.5–29	27–31	28.5–32.5	30–34.5	32–36.5	34.5–39.5	38–43
"C1"	16–18.5	17.5–20	19–22	20.5–24	22.5–26	24.5–28.5	27–31	28–32
"C"	14–16	15.5–18	17–20	18.5–21.5	20–23.5	22–25.5	24.5–28.5	25–29
Arm length/Longitud del brazo/Longueur de bras: (C-G) Regular: <48 cm Long: >48 cm								

Extra Wide circaid profile measuring chart

circum. (cm)	2	3	4	5	6	7	8	9
"G" extra-wide	N/A	32.5–36.5	34.5–39.5	36.5–41.5	39–44	41–46.5	44.5–50.5	48–54
"E" extra-wide	N/A	28–31.5	29.5–34	31–35.5	33–37.5	35–40	38–43	41–46
"C1"	N/A	17.5–20	19–22	20.5–24	22.5–26	24.5–28.5	27–31	28–32
"C"	N/A	15.5–18	17–20	18.5–21.5	20–23.5	22–25.5	24.5–28.5	25–29
Arm length/Longitud del brazo/Longueur de bras: (C-G) Regular: <48 cm Long: >48 cm								

